# MV-R

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
# Translation of Original Operating Instructions MV-R 100

  
EN

# Contents

<b>General information .....</b>	<b>4</b>
Limitation of liability .....	4
Responsibilities of the plant operator .....	4
Documentation .....	5
Content and structure .....	5
Marking system for integrated text and references .....	5
Formatting and symbols .....	6
Service address.....	6
<b>Safety.....</b>	<b>7</b>
Requirements of personnel .....	7
<b>Transport, packaging and storage.....</b>	<b>11</b>
Transport .....	11
Packaging.....	12
Storage.....	12
<b>Technical description .....</b>	<b>13</b>
Function .....	13
Conformity .....	14
Intended use.....	15
Possible misuse.....	15
<b>Technical data.....</b>	<b>16</b>
<b>Installation.....</b>	<b>17</b>
Positioning .....	17
Installation conditions .....	17
Electrical connection.....	18
Building installation.....	18
Rating plate .....	19
<b>Product description.....</b>	<b>20</b>
Main components of the product .....	20
Basic device .....	20
Back of basic device.....	21
Explanation of symbols.....	22
Interaction symbols.....	23
Unit symbols .....	24
Error symbols .....	25
Commissioning .....	25
Insertion of the door insulation .....	26
Explanation of main screens .....	27
System test.....	27
Error during system test .....	27
Main menu.....	28
Favourites.....	29
Overview of program options.....	29
Programs .....	30
New program .....	30
Service programs .....	32
My programs.....	32
Settings .....	33
Service .....	34
Generate QR code .....	34
Data transfer.....	35
Sintering aids.....	35
Preparations for sintering of zirconium oxide .....	35
Data transfer via USB stick.....	36
Restarting the device.....	37

On/off button on the touch display.....	38
<b>Maintenance and servicing.....</b>	<b>39</b>
Service programs .....	39
<b>Faults and error messages.....</b>	<b>40</b>
Safety .....	40
Overview of error codes .....	41
<b>Decommissioning.....</b>	<b>42</b>
<b>Disposal.....</b>	<b>42</b>
Safety .....	42
Disposal.....	42

## General information

### Limitation of liability

The content of these operating instructions has been drawn up in accordance with the applicable legislation and standards.

The device has been developed in accordance with the state of the art.

#### NOTE

**The manufacturer cannot accept any liability for damage resulting from:**

- Disregard/failure to observe the operating instructions
- Intentional misuse
- Use for a purpose other than that intended
- Operation by untrained personnel
- Operation by non-professionals (carrying out maintenance work, etc.)
- Technical modifications to the unit not agreed with the manufacturer
- Use of spare parts not approved by the manufacturer

### Responsibilities of the plant operator

The device is used in the commercial sector. The plant operator responsible for the device is therefore under the legal obligations of industrial health and safety. In addition to the safety information in these operating instructions, the safety, accident prevention and environmental protection regulations applicable to the area of application of the device must be observed.

#### The following applies in particular:

- The plant operator must familiarise himself with the applicable industrial health and safety provisions.
- The plant operator must ensure that all employees who use the device have read and understood these operating instructions.
- In addition, he must provide regular training for personnel and information about the risks that arise when handling the device.
- The plant operator must provide personnel with the necessary personal protective equipment.
- The plant operator must arrange for regular inspection of all safety equipment to ensure that it is complete and functioning properly.

## Documentation

### Content and structure

These operating instructions are an integral part of the device. They contain instructions and information concerning safe handling of the device and must be available to all users throughout the service life of the device. These operating instructions are aimed at trained operating personnel.

### Marking system for integrated text and references

The following types of notices are used:

#### **DANGER**

An imminent risk that could lead to serious physical injury or death.

#### **WARNING**

A potentially dangerous situation that could lead to serious physical injury or death.

#### **CAUTION**

A potentially dangerous situation that could lead to minor physical injury.

#### **NOTE**

A potentially harmful situation in which the product or an item in the immediate vicinity could be damaged.

#### **NOTE**

Note/tip for easier operation.

## Formatting and symbols

- △ *indicates a general piece of safety information*
- ☑ indicates that a requirement must be fulfilled
- 1. indicates steps in a series of actions
- ↻ indicates the result of an action
- indicates a list
- S1** indicates a button

## Service address

Friedrich-List-Straße 8  
D-76297 Stutensee-Blankenloch  
Tel.: +49 (0) 7244 70871-0  
[www.mihm-vogt.de](http://www.mihm-vogt.de)

## Safety

The **MV-R 100** is a high-temperature oven for commercial use in dental laboratories and may only be used to sinter sinterable ceramics.

### Requirements of personnel

Trained technical personnel who are familiar with handling the device, are able to carry out the work for which they are appointed and independently identify and avoid potential dangers on the basis of their technical training, knowledge, experience and familiarity with the relevant provisions.

#### NOTE

Personal protective equipment must be worn for all work on the **MV-R 100** in order to avoid accidents and damage to health.

#### WARNING

##### Electrical energy!

Risk of death from electric shock.

- Do not touch live cables and components with wet hands.
- Observe the accident prevention regulations for working with electric current.
- Before carrying out any installation, maintenance, cleaning or repair work, disconnect the power supply to the **MV-R 100** and secure it against being switched back on.

#### WARNING

##### Risk of ignition!

Use of flammable and explosive materials in the area of the oven.

- Do not operate the **MV-R 100** in the vicinity of flammable materials.
- Do not install the **MV-R 100** on flammable surfaces.

## WARNING

### **Risk of burning on hot surfaces!**

Surfaces become hot during operation of the *MV-R 100* and burns can occur if they are touched.

- Do not touch the housing or oven door during operation.
- Do not reach into the heating chamber. It can still retain high residual heat as a result of a preceding heating process.
- Allow the *MV-R 100* to cool down before carrying out maintenance, cleaning or repair work.
- Wear heat-resistant safety gloves if it is necessary to carry out work on hot components.
- Use suitable, sufficiently long tongs to insert and remove items during sintering.

## CAUTION

### **Incorrect operation!**

No liability can be accepted for damage that may be caused by misuse, incorrect operation, incorrect connections or improper maintenance/repair work carried out by untrained personnel. All warranty claims are also invalidated in such cases.

The unit must not be used if it or the mains cable becomes damaged and no longer functions correctly. In this case, contact the manufacturer immediately.

For your own safety and to increase the service life of your device, use only spare replacement parts.

To ensure safe operation of the *MV-R 100*, regional regulations (e.g. accident prevention regulations) apply in addition to these operating instructions. The former must be made available by the plant operator. The safety notices on the *MV-R 100* must be kept in a legible condition.

## CAUTION

### Materials and substances that are damaging to health!

Risk of contamination.

The insulation material in the combustion chamber is made of ceramic mineral fibres (INDEX no. 650-017-00-08). These are classified as carcinogenic (Category 1B (Annex VI, EC 1272/2008)). During work on the combustion chamber (e.g. replacement of the door insulation), fibre dust can be released.

The dust may have a carcinogenic effect if inhaled and cause irritation to the eyes, skin and respiratory organs.

- Wear personal protective equipment (safety glasses, protective gloves, respiratory protection of class FFP2, etc.) during all work on the combustion chamber.
- Use vacuum cleaners only with fine particle filters (HEPA dust class H).

In the area of the oven door, abrasion of insulation material can come about, which is visible as white dust. You should:

- wipe away these deposits with a damp cloth or
- clean them up with a vacuum cleaner with a fine particle filter (HEPA dust class H).

You must not blow away the dust deposits with compressed air.

First aid

- In the event of inhalation: If there is irritation, move to a dust-free area, drink water and clear your nose. If the symptoms persist, consult a doctor.
- In the event of skin contact: If there is skin contact, wash the affected area immediately with plenty of water and soap.
- In the event of eye contact: Immediately wash out your eyes carefully and thoroughly with an eye bath or water. Do not rub your eyes!

## NOTE

Always dispose of insulation material in a sealed container (see the "Disposal" section).

## WARNING

### Risk of burns due to objects falling out.

There is an acute risk of burns due to objects falling out. Ensure that the table and floor surfaces are fire-resistant.

1. For safety reasons, the oven must only ever be operated in an upright position. Only then is it possible for the operator to avoid objects that fall out quickly.
2. The plant operator at the operating location is responsible for providing appropriate safety training.
3. If there is a break in the sintering process because of failure to observe the recommended drying procedure, there is a risk of hot objects falling out as the oven door is opened.
4. Loading and unloading of the oven must be carried out exclusively with the aid of the tray fork. It is essential for safety reasons that a trained operator is present during the loading and unloading phase.

## NOTE

These operating instructions must be read and understood by each user before working on and with the device.

The operating instructions must be kept for the specified service life of the **MV-R 100**.

## Transport, packaging and storage

### Transport

#### WARNING

##### Injuries caused by *sintering oven falling!*

Slipping/falling when lifting and carrying the **MV-R 100** can lead to serious injuries.

- Carry/hold the **MV-R 100** only by the bottom edge of the housing (base).
- The **MV-R 100** must always be carried by at least four people.
- Mechanical strap eyelets are provided on the bottom of the device so that it can be carried more easily.

#### CAUTION

##### Risk of injury due to oven weight!

Physical overload/back problems due to heavy intrinsic weight.

- The **MV-R 100** must be carried/moved by at least four people.
- Mechanical strap eyelets are provided on the bottom of the device so that it can be carried more easily.

#### NOTE

##### Transport damage!

In order to avoid personal injury and transport damage:

- Transport the device in upright position only.
- Do not stack several devices on top of one another.
- Do not place any other objects on the device.
- Transport must be completed with as little shock and vibration as possible in order to avoid damage to the device.
- Ensure that the device is secured during transport to prevent it from slipping or falling over.
- If necessary, gloves (to prevent slipping and cuts) and safety shoes (to protect the feet in case the device is dropped) must be worn when carrying the device.
- The goods should be checked for damage and loss immediately on receipt and any damage confirmed on the delivery note so that a claim can be made against the freight forwarder. **Mihm-Vogt GmbH & Co.KG** cannot accept liability for damage and losses that are only identified at a later date.

## Packaging

### NOTE

The packaging protects the *MV-R 100* against transport damage, corrosion and other damage. Therefore remove it only just before commissioning and store it in a dry place for later use.

## Storage

### NOTE

#### Temperature damage!

To avoid temperature damage:

- Store the device only between +5°C and +40°C.
- Always store the device in a dry, dust-free location.
- Avoid direct exposure to sunlight.
- Avoid mechanical shocks.

## Technical description

### Function

The **MV-R 100** is used to process sinterable ceramics.

The product to be sintered is placed in the sintering tray and positioned on the support plate. Once the heating parameters have been entered and the start symbol pressed, the electrically operated oven door closes and the heating process begins.

When the heating program has been completed and the **MV-R 100** has cooled down, the oven door opens and the finished product can be removed.

### Heating chamber

The product is sintered in the heating chamber. This consists of two different layers of ceramic insulation and is operated with  $\text{MoSi}_2$  heating elements. The outer insulation layer is designed for temperatures up to  $1200^\circ\text{C}$ , the inner layer for temperatures up to  $1700^\circ\text{C}$ .

### Oven door

The oven door has a ceramic door filling.

A slip clutch in the drive mechanism prevents excessive contact pressure between the oven door and heating chamber.

### Oven housing

The oven housing is made of steel plate and is plastic-coated inside and out.

### Program controller

The program controller has a timer function for day of the week and time.

Operating parameters and heating programs are stored in a non-volatile memory and are retained even if there is a power cut.

The target temperature set is maintained with an accuracy of  $\pm 1^\circ\text{C}$ .

A temperature sensor integrated into the heating chamber records the chamber temperature near the product.

A thermocouple breaking device prevents overheating of the **MV-R 100** as a result of a defective temperature sensor.

## Conformity

### EC declaration (Translation of the original declaration of conformity)

The manufacturer/distributor

MIHM-VOGT GmbH & Co. KG  
Friedrich-List-Str. 8  
76297 Stutensee  
Tel.: +49 (0) 72 44/7 08 71-0  
Fax: +49 (0) 72 44/7 08 71-20  
Email: info@mihm-vogt.de

hereby declares that the following product

Product designation: Sinterofen  
Model name: MV-R 100  
Year of construction: 2020

**Description:**

The MV-R 100 is a high-temperature oven for commercial use in dental laboratories and may only be used to sinter sinterable ceramics.

fulfills all the relevant provisions of the legislation (in the following) - including the changes which applied at the time of the declaration. This declaration of conformity is issued under the sole responsibility of the manufacturer. This declaration only refers to the machine in the condition in which it was placed on the market; it does not apply to parts installed by the user and/or modifications made at a later time.

The following legislation have been applied:

Machinery Directive 2006/42/EC  
EMC Directive 2014/30/EU  
RoHS Directive 2011/65/EU

The following harmonised standards have been applied:

EN 61010-1:2010  
EN 61010-2-010:2014  
EN 61326-1:2013  
EN ISO 12100:2010

Name and address of the person authorized to compile the technical file:  
Tobias Gillen

Place: Stutensee  
Dat: 29.10.2020

## Intended use

The **MV-R 100** is a high-temperature oven for commercial use in dental laboratories and may only be used to sinter sinterable ceramics.

### NOTE

In the event of misuse, incorrect operation, incorrect connections or improper maintenance/repair work carried out by untrained personnel, safety of personnel can no longer be guaranteed and no liability can be accepted for any resulting damage. All warranty cover is also excluded in such cases.

Use of spare and wear parts that were not purchased from the manufacturer invalidates the certification and warranty for the **MV-R 100**.

## Possible misuse

- Deployment of personnel who have not been trained and are not adequately qualified.
- Use of products that are not approved by the manufacturer.
- Use of spare parts that are not approved by the manufacturer.
- Use that is not in accordance with the Declaration of Conformity.
- Technical modifications and changes to the device that have not been approved by the manufacturer.

### NOTE

Only sintering aids/accessories/wear parts and spare parts approved by Mihm-Vogt may be used.

## Technical data

<b>MV-R 100</b>		
<b>General information</b>	MV-R 13A	MV-R 15A
Dimensions [mm] (W x D x H)	390 x 540 x 780 mm	
Combustion chamber volume	2 trays Ø 100/30mm	
Max. temperature [°C]	1650°C	
Weight [kg]	65 kg	
Minimum gap around the device	100 mm	
<b>Electrical connection values</b>		
Power supply	220 - 230 V	200 - 240 V
Frequency	50 Hz	50/60 Hz
Max. power consumption	3500 W	
Overvoltage category	Overvoltage category II	
<b>Fuse protection</b>		
In device	16 AT	
Provided by the customer	Connection to a separate electrical circuit with a 16 A type K, Z fuse (other types of fuse in accordance with the regulations in the country of use)	
Protection type	IP 20 (protection against penetration of foreign bodies but not of water)	
<b>Operating conditions</b>		
Installation area	Indoors only (in dry rooms)	
Temperature range	+5 / +40°C	
Relative humidity	Up to 31°C: 80%	
Maximum humidity	Up to 40°C: 50% no condensation	
Altitude	Max. 3000 m	
Degree of contamination	2	

# Installation

## Positioning

The *MV-R 100* is designed as a table-top device. For stable positioning, a level surface of at least 50 cm x 60 cm with a load-bearing capacity of 80 kg is recommended.

### Installation conditions

- ▶ Always set up the *MV-R 100* in a dry room that is as dust-free as possible and ensure that no liquids can come into contact with the device.
- ▶ There must be no easily ignitable or flammable gases or liquids stored in the installation rooms.
- ▶ Do not put any flammable or ignitable objects in the vicinity of the *MV-R 100*.
- ▶ Keep a minimum gap of 100 mm around the *MV-R 100* to ensure adequate cooling.

### CAUTION

#### Risk of ignition!

Ignition of the floor covering due to product tray breaking.

- ▶ Ensure that the floor covering does not have a flammable surface.

### CAUTION

#### Tipping loads!

Inadequate load-bearing capacity of the installation surface.

Ensure that the surface has adequate load-bearing capacity when setting up the *MV-R 100*.

### CAUTION

#### Risk of injury due to oven weight!

Physical overload/back problems due to heavy intrinsic weight.

- ▶ The *MV-R 100* must be carried/moved by at least four people.
- ▶ Mechanical strap eyelets are provided on the bottom of the device so that it can be carried more easily.

### CAUTION

#### Risk of overheating!

Overheating due to blocked air vents.

- ▶ Ensure that the ventilation slots remain free on all sides.

1. Make sure the installation surface is horizontal.
  2. Place the **MV-R 100** on the installation surface.
- ⚠ *Ensure that the surface is non-slip.*
- ⚠ *Lift and carry the sintering oven only by the bottom of the device.*

## Electrical connection

### Building installation

- The **MV-R 100** requires its own electrical circuit.
- The electrical circuit in the building must have a fuse device of minimum 16 A type K, Z (other fuse types in accordance with regulation in the country of use).
- An additional circuit breaker (with a tripping current of 30 mA) must also be installed.
- For safe electrical operation, the **MV-R 100** requires a protective conductor connected to the socket.
- When selecting the installation location, it must be noted that the mains cable supplied is 2.0 m long and cable extensions are not permitted. The power supply must be in the nominal voltage range of 200 - 240 volts.


### DANGER

#### Electrical energy!

Risk of death from electric shock.

- Do not touch live cables and components with wet hands.
- Observe the accident prevention regulations for working with electric current.
- Connect the device only to a power supply that matches the information on the rating plate.
- Removable mains cables must not be replaced with mains cables of inadequate size.

### Rating plate


- |   |                          |   |                            |
|---|--------------------------|---|----------------------------|
| 1 | Machine type/designation | 6 | Mihm-Vogt reference number |
| 2 | Operating voltage        | 7 | Serial number              |
| 3 | Mains frequency          | 8 | Year of manufacture        |
| 4 | Output                   | 9 | Manufacturer's details     |
| 5 | CE mark                  |   |                            |

## Product description


### Main components of the product

#### Basic device

- 1 Basic device
- 2 Control panel
- 3 Display on button
- 4 Status display
- 5 Door insulation/base frame
- 6 Door insulation overlay
- 7 Sinter trays (not shown)

## Back of basic device

- 1 Rating plate
- 2 Main switch on/off
- 3 Mains connection
- 4 USB port  
(power supply: 5 V, rated current: 500 mA)

**Note:**

The USB port is used exclusively for data transmission, back-ups and software updates on USB stick.

No other devices – e.g. power banks – may be connected.

- 5 LAN connection
- 6 Main fuse F1
- 7 Main fuse F2
- 8 Control thermocouple (not active)

## Explanation of symbols

The function is active when the symbols turn yellow.

### Operating symbols

	Main control system menu
	Continue to next possible step
	Back to last step
	Start selected program
	Stop running program
	Open/close oven door
	Activate/deactivate a timer function
	Activate/deactivate a pre-heating function
	Activate/deactivate a pre-drying function
	SpeedPlus program pre-selection
	SpeedMotion program pre-selection
	Selection of a saved program on the "Set programs" or "My programs" lists and display on "My list"

### Interaction symbols

Confirmation of an interaction that appears as a pop-up

Cancellation of an interaction that appears as a pop-up

Info: Error present

12H display AM/PM

24H display

Keyboard input


EU date display

USA date display

Keyboard/input field

**Unit symbols**

Increase rate in °C/min

Set temperature in °C

Holding time in min

Target time in hh:mm

Actual temperature in °C

**Results symbols**

Display when process has finished without errors

Display when process has finished with errors

**Warning symbols**

Warning of heat when opening the door during the process

### Error symbols

Informative error (4-digit error code) process continues  
Error codes can be displayed by tapping

Less serious or serious error (2 or 3-digit code)  
Error codes can be displayed by tapping

Less serious error (3-digit error code)  
Process continues/work must be checked

Serious error (2-digit error code)  
Automatic cancellation of process/work is very likely to be unusable

### Commissioning

#### NOTE

In temperatures below 15°C (e.g. after transport), allow the machine to stand at room temperature for approx. 60 minutes before commissioning.

1. Connect the power supply cable to the mains.
2. Switch the oven on at the main switch.
- ↻ The door of the device opens and the message "Insert door insulation" appears.
3. Insert the door insulation and confirm the message with a tick.

## Insertion of the door insulation

- 1 Door insulation base support mv-r, 100/70
- 2 Ceramic protection tube, 6/3x90
- 3 Door insulation overlay mv-r, 100/40
- 4 Support fleece, 100; for sinter tray

- The device then carries out a self-test and starts in commissioning mode.
- 4. Follow the commissioning menu guidance
- On completion of commissioning, the device is on the main menu with the door open.

## Explanation of main screens

### System test

The system test is carried out when the device is started. Wait for this phase to finish and do not switch the device off.

Successfully completed tests are marked with a green tick. The green bar indicates the progress of a test.

### Error during system test

If there is an error, the corresponding error code is displayed. A red cross appears after the defective function. Follow the instructions on the screen.

#### NOTE

If functions are defective, repeat the self-test. If the problem persists, please contact customer service.

## Main menu

The following items appear on the “Main menu” screen and can be selected:

- Favourites
- Program options
- Settings
- Service
- Generate QR code
- Open/close door
- Data transfer

The “Service” section is only for the service technician.

## Favourites

You can switch to the “Favourites” screen from the “Main menu” screen. This shows an overview of your favourite programs.

You can remove these programs from the list. To do so, tap on the symbol next to the program in question.

On the “My programs” screen, you can add existing programs to your favourites or remove them. You can select your favourite programs here and start them directly.

## Overview of program options

You can switch to the “Programs” screen from the “Main menu” screen. The following items appear on the menu:

- Programs
- New program
- Service program
- Favourites

## Programs

Under the Programs menu item, you will find an overview of all program lists saved to the system. These include:

- My programs
- SpeedPlus sintering programs 
- SpeedMotion sintering programs 

## New program

You have 3 options for creating new programs:

- ▶ Standard sintering programs (< 120°C/min)
  - Allows input of heating rates from 1°C/min to 120°C/min
  - Up to 10 steps freely programmable (S1 to S10)
  - Linear cooling (lift opens 30% from a temperature of 1100°C)
- ▶ SpeedPlus sintering programs 
  - Allows input of heating rates from 1°C/min to 200°C/min in pre-defined temperature ranges which are default maximum values for steps S1 to S3.
  - Up to 7 steps freely programmable (S4 to S10)
  - Linear cooling (lift opens 100% from a temperature of 1650°C)
- ▶ SpeedMotion sintering programs 
  - Pre-set max. possible heating rates
  - Consists of just 1 step, only sintering temperature and holding time can be programmed.
  - At the end of the holding time, the process is complete and the lift opens 100%.

You also have the opportunity to pre-set the following program options before each sintering program:

- ▶ Pre-drying option 
  - Each time this symbol is pressed, a pre-drying time of 15 mins. is set. Pressing the symbol multiple times allows periods of 30/45/60 mins. to be set. (Visible from the number of heating waves that turn yellow)
  - When the sintering program is started with this option selected, the lift stays 30% open and the restoration can dry for the set time at approx. 80°C. At the end of this time, the lift closes and the programmed sintering program starts directly after.
  
- ▶ Pre-heating option 
  - This option is recommended for fast sintering programs in particular.
  - Immediately on selection of this option, the lift and the oven heats to 250°C. This option ends after 30 mins. without any interaction with the user. Pre-heating can be cancelled actively by the user to load the oven and then start the programmed sintering program manually.

You can switch to the "New program" screen from the "Programs" and "My programs" screen. New programs can be created here.

1. To create a new program, tap on the first field – Increase rate – on the next free row.
  - ↻ The border of the field turns orange.
2. Enter the increase rate in °C/min using the number pad.
3. Tap on the Enter button to complete the input.
  - ↻ The orange border jumps to the next field: "Temperature".
4. Enter the temperature in °C using the number pad.
5. Tap on the Enter button to complete the input.
  - ↻ The orange border jumps to the next field: "Duration of first process stage".

6. Enter the duration of the first process stage in minutes using the number pad.
7. Tap on the Enter button to complete the input.
- ↻ The orange border jumps to the first field on the next row. Creation of the program is complete.
8. After creating the program, tap on the "Next" button.

## Service programs

You can find further information in the "Maintenance and servicing" section under "Service programs".

## My programs

You can switch to the "My programs" screen from the "Programs" screen. This shows an overview of recent programs.

You can add these programs to your favourites or remove them. To do so, tap on the symbol next to the program in question. The symbol then changes colour from white to dark grey or dark grey to white.

A grey symbol means that the program is not available as a favourite.

A yellow symbol means that the program is available as a favourite.

## Settings

### Language

Tap here to set the language for the menu guidance.

### Temperature unit

Tap here to set the temperature unit.

### Date & time

Tap here to set the date & time.

### Device details

Tap here to obtain information about the model/series, serial number, hardware/software version and operating hours.

### Factory settings

Tap here to reset the device to its factory settings.

## Service

### Device details

Tap here to obtain information about the model/series, serial number, hardware/software version and operating hours.

### System test

Tap here to start the system test manually.

### Error

Tap here to display and, if necessary, clear the error list.

### Service programs

Tap here to run the service programs.

## Generate QR code

Tap here to generate a QR code on the device display containing the following information:

- Model/series
- Serial number
- Hardware/software version
- Operating hours
- All undeleted, saved error codes

## Data transfer

Tap here to activate data transfer by USB stick.

## Sintering aid

### NOTE

Only use sintering aids approved by Mihm-Vogt.  
Consult the information leaflet for the sintering aid in question for application information.

### CAUTION

#### Risk of crushing of limbs!

The oven door closes automatically.

- Press the **Start** button only when the filled sintering tray has been placed in position.
- Ensure that no one reaches in between the oven door and heating chamber as the oven door is closing.

## Preparations for sintering of zirconium oxide

Recommendation for filling the sintering tray:

**Note: The SiC sinter tray can only be used up to a maximum temperature of 1550°C!**

1. Fill the sintering tray with a whole bottle of sintering granules.
2. Place the parts to be sintered with the chewing surface facing downwards on the sintering granules. Ensure that the parts to be sintered are not touching one another or the sintering tray.
3. Insert the base frame (1) into the door, put the fixing pin for the door insulation (2) in the centre hole provided for it in the base frame (1), place the support brick (3) on the fixing pin (2) and finally place the fleece (4) in the centre of the support brick (3) on the guide pin.
4. Place the sintering tray on top.
5. You can load and unload the tray with the removal tongs.

## Data transfer via USB stick

To transfer data to the device via USB stick, a USB stick 2.0 without password protection must be provided.

1. Connect the USB stick to the back of the device.
2. The USB stick is automatically recognised as a storage medium and you have the following options:


- “Firmware update”  
Tap here to update the device firmware.
  - “Update configuration”  
Tap here to update the device configuration.
  - “Load my programs”  
Tap here to transfer the programs on the stick to the device.
  - “Back up my programs”  
Tap here to back up the programs you have created.
  - “Back up data”  
Tap here to create a complete back-up of the device.
  - “QA document”  
Tap here to create a QA document for the last 10 completed processes. You will then find the document on the USB stick.
3. Following successful transfer of the data, remove the USB stick and carry out a restart.

## Restarting the device

To carry out a restart, switch the device off at the back and wait approx. 30 seconds until the green light stops flashing and goes out.

Switch the device on again at the main switch.

## On/off button on the touch display

To protect the touch display, switch the display off via the on/off button. You must press and hold the button for at least 1 sec. to do this. The device is ready for the next sintering process and does not have to be restarted again when the display is switched on.

## Maintenance and servicing

Clean the housing of the *MV-R 100* now and again with a damp cloth.

### NOTE

#### Damage to the heating!

- Ensure that the heating chamber does not become dirty.  
The heating could be damaged.

### NOTE

#### Impairment of service life due to colouring liquids!

- During the sintering process, colouring liquids can shorten the life cycle of the heating elements significantly.

## Service programs

A cleaning cycle should be run in accordance with the frequency of use. This is used to remove contamination caused by liquids and other dirt that is deposited in the insulation.

A cleaning regeneration cycle should also be run in accordance with the frequency of use, which is necessary to regenerate the heating elements.

- Temperature check (only in conjunction with test kit)
- Cleaning of heating chamber
- Regeneration of heating elements

## Faults and error messages

### Safety

#### DANGER

##### Electrical energy!

Risk of death from electric shock.

- Work on electrical systems may be performed by qualified electricians only.
- Before carrying out any installation, maintenance, cleaning or repair work, disconnect the power supply to the **MV-R 100** and secure it against being switched back on.
- Do not touch live cables and components with wet hands.
- Observe the accident prevention regulations for working with electric current.

#### WARNING

##### Hot surfaces!

Serious burns to the limbs.

- Do not touch the housing or oven door during operation.
- Allow the **MV-R 100** to cool down completely before carrying out maintenance, cleaning or repair work.
- Wear heat-resistant, heat-insulated safety gloves if it is necessary to carry out work on hot components.

#### NOTE

##### Material damage as a result of defective repairs to electrical cables!

Malfunctions and defective electrical components possible.

- Do not repair defective cables and plugs.

### Overview of error codes

Error	Category	Cause	Measures
F10	Control system	GS1 defective	Service/dealer
F1000	Control system	SD card defective	Service/dealer
F1010	Control system	Battery defective	Service/dealer
F20	Mains voltage	Mains supply interrupted	Check mains supply
F2010	Mains voltage	Mains supply fault	Check mains supply
F30	Heating circuit	GS1 defective (no current flow in power unit)	Service/dealer
F32	Heating circuit	GS1 defective (power unit error message)	Service/dealer
F40	Thermal sensors	Temperature sensor S short circuit	Service/dealer
F41	Thermal sensors	Temperature sensor S sensor break	Service/dealer
F42	Thermal sensors	Thermal sensor on circuit board defective	Service/dealer
F400	Thermal sensors	Temperature sensor K short circuit	Service/dealer
F401	Thermal sensors	Temperature sensor S sensor break	Service/dealer
F50	Cooling	Cooling system defective	Service/dealer
F500	Cooling	Fan defective	Service/dealer
F60	Drive	Door drive, door cannot be positioned	Service/dealer
F61	Drive	Door drive, lower limit switch defective	Service/dealer
F62	Drive	Door drive, upper limit switch defective	Service/dealer
F70	Process temperature	Temperature error	Usability of the sintered product, check power supply
F71	Process temperature	Maximum temperature above 1680°C emergency shut-off	Service/dealer
F700	Process temperature	Process temperature outside tolerance	Usability of the sintered product, check power supply
F701	Process temperature	Final temperature not reached	Usability of the sintered product

## Decommissioning

Decommissioning can be carried out for two reasons:

- With the aim of setting up again at a different location.
- With the aim of final disposal.

If the **MV-R 100** is to be set up at a different location, decommissioning must be prepared well. All components and fittings must be disassembled carefully, labelled and, if necessary, packed for transport. This ensures when it is set up again that all parts are correctly allocated and can be re-installed in the correct place.

1. Switch the **MV-R 100** off.
2. Disconnect the **MV-R 100** from the power supply.
3. Disconnect all connections from the **MV-R 100**.

## Disposal

### Safety

#### WARNING

##### Release of harmful substances!

Harmful substances can be inhaled when the insulation materials are handled.

- Personal protective equipment (respiratory protection) must be worn during disposal.

#### WARNING

##### Contamination of the environment and groundwater due to improper disposal!

- The statutory regulations and guidelines in the country of operation must be observed when disposing of parts of the unit and operating materials.

### Disposal

1. Sort the component parts of the **MV-R 100** into recyclable materials, hazardous substances and operating materials.
2. Dispose of the components of the **MV-R 100** or take them to be recycled.

# CONSISTENT. RELIABLE. PRECISE.

**MIHM-VOGT GmbH & Co. KG**  
Friedrich-List-Straße 8  
76297 Stutensee-Blankenloch  
Germany

Web [www.mihm-vogt.de](http://www.mihm-vogt.de)  
Mail [info@mihm-vogt.de](mailto:info@mihm-vogt.de)  
Phone +49 7244 70871-0  
Fax +49 7244 70871-20