


D-one
What you see is what you get!

D - ONE



Technical parameters

Name:D-one milling machine

Weight:80kg

Max power:900W

Number of axis:4axis

Processing mode:wet

Max feed rate:6000mm/min

Tool changing mode:automatic

Processing materials:PMMA and Glass ceramic

Dimension (L*W*H):700*490*610(mm)

Input voltage:AC220V

Spindle power:250W

Travel range: 360° (positioning processing at any angle)

Max speed:60,000rpm

Tool magazine capacity:3

Spindle cooling:air cooling

Processing time:Glass ceramic:25mins



Professional Service Team



Online Support



Seasonal On-site Technical Support



Add: 2nd Floor, 1st building, No.88, North Zhuguang Rd, West Area, Honghualing Industrial Park, Nanshan District, Shenzhen, China

Tel: +86-755 8600 1804
E-mail: xtcera@xianton.com

Fax: +86-755 8600 1486
Web: en.xtcera.com