





MACHINE SETUP



MACHINE SETUP DC5™ / DC7™

So that materials such as PMMA, composite, PEEK or similar materials can be handled by the filter system of the DC5™ and DC7™, a new filter box has been developed.

Certain types of debris can otherwise block the valves, filters or supply lines during wet processing. The debris can now be separated out in advance with the aid of the new filter box.




DC5™ / DC7™ filter box:

- ✓ New filter membrane separates debris quickly and simply
- ✓ Simple installation and removal
- ✓ Clever and inexpensive
- ✓ Usable in all standard machines from the DCS series
- ✓ Also suitable for the separation of CoCr or titanium debris
- ✓ No additional replaceable filter
- ✓ Can be cleaned quickly
- ✓ Prolongs the life of the cooling fluid and filter
- ✓ Made in Germany

Price per DC5™ / DC7™ filter box::

119,50 €*



DC5™ / DC7™ filter mat set:

A filter mat set is offered for the grinding of ceramic materials. The filter mats can be placed in the debris drawers. The grinding dust is trapped in the mats and can be quickly disposed of.

Price per DC5™ / DC7™ filter mat set:
(contains 2 filter mats)

16,40 €*

* Prices exclude tax and freight charges. Offers subject to change without notice.

CAD/CAM Hotline: + 49 (0) 731 14661111

www.dental-concept-systems.com



DC5™ DC7™



All parts from our own production – made in Germany

I am interested in the following product:

☐

DC5™ / DC7™ filter box

☐

DC5™ / DC7™ filter mat set

My/our contact details:

Company:

Contact:

DC5 serial number (*label on the back of the machine*):

Street/number:

Post code/place:

Phone:

Email:

Notes/remarks:

Please email to: service@dental-concept-systems.com

Distributed by:



4220 Steve Reynolds Blvd. Ste 24
Norcross, GA 30093

GA Office: (770) 921-2131
TX Office: (972) 600-9187
sales@zublerusa.com

