


AIR HOSE RETURN DEVICE with 6/3 mm hose, air blower with air adjustment valve, connector hose and air socket

Item No. | 20148/20238


Space-saving, gentle tube routing with definitive stop (no pull) during use. Compact unit, for installation anywhere.

TECHNICAL DATA	20148	/	20238
Height: 90 mm, Depth: 30 mm			
Length:	420 mm	/	470 mm
Tube extension:	550 mm	/	690 mm
Connector:	8/6 mm	/	8/6 mm
Discharge:	6/3mm	/	6/3mm
Weight:	700g	/	750g



■ | Air hose return device, with 6/3 mm hose 420 mm

■ | Air hose return device, with 6/3 mm hose 470 mm



■ | Holder straight for fitted air blower (drill hole \varnothing 33 mm)

■ | Holder slanted 30° for fitted air blower (drill hole \varnothing 33 mm)

■ | Holder straight for air blower, for sub-structure assembly

■ | Holder slanted 30° for air blower, for sub-structure assembly



■ | Air blower with connector for 4/3 mm hose

■ | Air blower with connector for 6/4 mm hose


■ | Air blower with connector for 8/6 mm hose

RETURN DEVICE for handpiece cable

Item No. | 20178/20239

A space saving gentle tube guidance with definitive stop (no pulling) during use. Compact unit, for installation anywhere.

TECHNICAL DATA	20178	/	20239
Height: 90 mm, Depth: 30 mm			
Length:	420 mm	/	470 mm
Tube extension:	550 mm	/	690 mm
Weight:	500g	/	550g



■ | Holder for handpiece \varnothing 28 mm (drill hole \varnothing 37 mm)

■ | Holder for handpiece \varnothing 30 mm (drill hole \varnothing 37 mm)

■ | Holder for handpiece \varnothing 32 mm (drill hole \varnothing 37 mm)