




VARIOstar T50

Tabletop handpiece system



Distributed by:

zubler USA

Zubler USA LLC
805 Marathon Pkwy. Ste 130
Lawrenceville, GA 30046

GA Office (770) 921-2131
TX Office (972) 600-9187
www.zublerusa.com





System-Variants:

804 - 3112 VARIOstar K50 Knee control unit 220 V - 240 V

804 - 3111 VARIOstar K50 Knee control unit 115 V

Contains:

804 - 3012 Knee control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3011 Knee control unit 115 V



804 - 3152 VARIOstar T50 Tabletop unit with Foot pedal 220 V - 240 V

804 - 3151 VARIOstar T50 Tabletop unit with Foot pedal 115 V

Contains:

804 - 3052 Tabletop control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3051 Tabletop control unit 115 V

804 - 1016 Foot pedal



804 - 3162 VARIOstar T50 Tabletop unit with Foot control lever left/right 220 V - 240 V

804 - 3161 VARIOstar T50 Tabletop unit with Foot control lever left/right 115 V

Contains:

804 - 3052 Tabletop control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3051 Tabletop control unit 115 V

804 - 1015 Foot control lever



804 - 3172 VARIOstar K50 Tabletop unit with Knee pedal 220 V - 240 V

804 - 3171 VARIOstar K50 Tabletop unit with Knee pedal 115 V

Contains:

804 - 3052 Tabletop control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3051 Tabletop control unit 115 V

804 - 1017 Knee Pedal (left/right mountable)



All systems include a tool for changing the collet chuck, a handpiece holder and operating instructions.

Technical data:

Length:	162 mm	Variants:	Knee control unit
Weight:	235 g		Tabletop unit with knee pedal
Rated voltage:	AC 100V/115V/230V/240V		Tabletop unit with foot pedal
Rated frequency:	~ 50/60Hz		Tabletop unit with foot control left/right
Speed:	1000-50000 rpm	Warranty:	Lifetime warranty on armature and winding
Chuck:	Ø 2.35 mm or optional Ø 3.00 mm		All other parts 1 Year
Torque:	7,8 Ncm		
Power consumption:	max. 270 W		
Protection class:	I		

Subject to technical changes.

Distributed by:



Zubler USA LLC
 805 Marathon Pkwy. Ste 130
 Lawrenceville, GA 30046

GA Office (770) 921-2131
 TX Office (972) 600-9187
www.zublerusa.com

