

## Articles in this section



# System requirements for the i500

1 month ago · Updated

[Follow](#)

## [1] i500 Intraoral Scanner System Requirements

The following are the system requirements in order to use the i500 intraoral scanner with the iScan software. We recommend checking that your system meets the minimum requirements for smooth operation.

### - Laptop

- CPU: Above Intel Core i7-8750H
- RAM: Above 16 GB
- Graphics Card: Above NVIDIA GeForce GTX 1060
- OS: Windows 10 64-bit
- Use computer and monitor certified IEC 60950, IEC 55032, IEC 55024

### - Desktop

- CPU: Above Intel Core i7-8700
- RAM: Above 16 GB
- Graphics Card: Above NVIDIA GeForce GTX 1060
- OS: Windows 10 64-bit
- Use computer and monitor certified IEC 60950, IEC 55032, IEC 55024

You can refer to the Medit i500 User Guide if you require additional information.

1. Run Medit iScan.
2. Click the **Menu** button at the top left corner.

