

A2 LED

1000 HOUR LED LIGHT GUARANTEE

New Arrival





CONTENTS

01



A2 Advantages

02



Compare

03



Specials





sales@zublerusa.com



ADVANTAGES

A2

LED Handpiece


Impervious LED Technology

Our SMD LED Technology allows the hand piece to be autoclaved more than 1200 times under 135°C and includes a 1000 Hour working life guarantee!





ADVANTAGES



Strong 22 Watt
High Torque Motor



SMD LED
Technology



Perfect Cooling
Effect



350,000 RPM



Japanese Ceramic
Covered Bearings




Quiet
Less Than 65dB

SMD LED Technology

≥1200 times autoclaved


≥1000 hours LED guarantee



sales@zublerusa.com

Integrated Solid Mold Head

Provides smooth spray without clogging
and no risk of breakage at the head and
neck!





ADVANTAGES

Technical Data:

Model Version: A2

Model TUP

Working Pressure: 0.28(Mpa)-0.3(Mpa)

Rotational Speed: 300(krpm)-350(krpm)

Water flow: 60(ml/min)-90 (ml/min)

Noise: ≤ 65 dB

Light: White SMD LED

Brightness: 16,000Lux-20,000Lux

Color Temperature: 5,500(K)-6,000(K)

Head Diameter 12.5mm

Height: 13mm

Description: Torque Head, push button, triple water spray, exposed cartridge.



A2

LED handpiece



Distributed by:



805 Marathon Parkway
Suite 130 Lawrenceville,
GA 30046

(770) 921-2131

16135 Preston Rd.
Suite 126
Dallas, TX 75248
(972) 600-9187

www.zublerusa.com

sales@zublerusa.com

