

Make good press results perfect press results!

This intelligently designed device from Zubler is the solution to surface anomalies in pressable ceramics!

- Immediately after removing a ring from the oven, place it into the **VARIO SC Surface Conditioner**.
- The **VARIO SC** adds the necessary force to compensate for the contractions that occur during the cooling process, thus avoiding undesirable malformation.
- Furthermore, the cooling fan and stability added by the **VARIO SC** reduces cooling time by as much as 50%!



Order these also...

CONCEPT VEST™

The special investment for lithium disilicate ceramics

The investment specifically developed for processing lithium disilicate glass ceramics guarantees ...

- Perfect surfaces
- A minimal reaction layer after pressing
- Detailed reproduction of the pressed objects
- Consistent expansion values through strictest quality controls



25 x 100 g
 50 x 100 g
 Mixing Liquid / 1 Liter

Item-No. 898/2066
Item-No. 898/2065
Item-No. 898/204

Zubler® Disposable Plungers

For pressing of all common dental ceramics
 and **lithium disilicate**.

Their excellent quality and precise dimensions guarantee
 process control and consistent results for any ring size.



ø 12 mm **Item-No. 20/0235**
ø 13 mm Item-No. 20/0238
 ø 16 mm **Item-No. 20/0237**

Zubler® Investment System for Lithium Disilicate Ceramics

Zubler Flex-ring formers create extrusions on the surface of rings ensuring faster heating and cooling. Sector lines make it easier to locate pressed objects for divesting. The Flexing is easy to clean and has a much longer service life than most other ring formers.



Flexring 100 g Set **Item-No. 20/0228**
 Flexring 200 g Set **Item-No. 20/0229**



Zubler USA
 805 Marathon Parkway, Ste 130
 Lawrenceville, GA 30046

Georgia Office: (770) 921-2131
 Texas Office: (972) 600-9187
www.zublerusa.com



Distributed By:

 **zublerUSA**