



# THE CLEAN SOLUTION

Catalog | [www.bdt-dental.de](http://www.bdt-dental.de)

# WHEN QUALITY HAS TRADITION AND OUTSTANDING SERVICE IS, WHAT MAKES US SUCCESSFUL

## How do you go from a small company in Offenburg to a globally successful manufacturer?

Ditmar Stemmer, Anton Prinzbach und Klaus Bruder founded BDT GmbH on the 1<sup>st</sup> January 1983, carrying on the dental technology division as a separate business when Bregler KG ended equipment manufacturing in 1982. Ditmar Stemmer continues as a Managing Director to this day following the retirement of the other founders.

The mission was clear: Improve work in dental laboratories by simplifying and optimizing processes and end products. In order to achieve this, we developed our own tool molds and had the necessary components manufactured by renowned, regional companies, which were then assembled in our facilities in Offenburg. This way, we could ensure the highest quality and efficiency for the production of larger quantities in a shorter period of time – in those days as well as today.

The first in-house product developments included the polymerization bath, the plaster separator with the 4-chamber clarification system and the air hose return device. The interest was great and the success confirmed our mission. Our products were in increasing demand, so that from now on we are regularly expanding the range of dental products. We have also been ISO-9001 certified since 2019.

Today, with Jessica Cooper, Ditmar and Daniel Stemmer as managing directors, we are a family business in the second generation, we sell our sophisticated dental products worldwide and are an international, successful manufacturer with the highest quality “Made in Germany”.

As a family business we stand for informative, fast and flexible interaction in all areas. With BDT at your side, you benefit from 40 years of experience and top quality.

*Ditmar Stemmer*

**Ditmar Stemmer**  
Managing Director

*Jessica Cooper*  
**Jessica Cooper**  
Managing Director

*Daniel Stemmer*

**Daniel Stemmer**  
Managing Director



# POLYMERIZATION TECHNOLOGY

Polymerization bath VO - 2	44
Polymerization bath VO - 4	45
Cooling baths	46
SR IVOCAP accessories	47
Wall holder	48
Stand holder	49







# POLYMER- IZATION TECHNOLOGY

Good teeth and strong denture are closely related to one's own quality of life. That is why we ensure precise and high-quality prostheses with our polymerisation technology. The focus is always on high wearing comfort and an accurate fit. In our range you will find everything you need to put a radiant smile on your patients – comfortable and radiant.

# POLYMERIZATION BATH

for 2 cuvettes

- ✓ Stainless steel bath
- ✓ Stepless continuous temperature regulation
- ✓ Indirect heating

Item no. 10002

IVO-2



Item no. 10004

IVO-2 special



Item no. 10005

IVO-2 manual



## PRODUCT INFORMATION

	IVO-2	IVO-2 special	IVO-2 manual
Dimensions (l x w x h)	350 x 470 x 300 mm	350 x 360 x 260 mm	350 x 360 x 260 mm
Electrical values	230 V / 50 Hz / 1,500 W	230 V / 50 Hz / 1,500 W	230 V / 50 Hz / 1,500 W
Compressed air	X	–	–
Timer	X	–	–
Water level control	X	X	–
Item no.	10002	10004	10005



# POLYMERIZATION BATH

for 4 cuvettes

- ✓ Stainless steel bath
- ✓ Stepless continuous temperature regulation
- ✓ Indirect heating

Item no. 10003

IVO-4



Item no. 10033

IVO-4 special



Item no. 10006

IVO-4 manual




## PRODUCT INFORMATION

	IVO-4	IVO-4 special	IVO-4 manual
Dimensions (l x w x h)	530 x 470 x 300 mm	530 x 360 x 260 mm	530 x 360 x 260 mm
Electrical values	230 V / 50 Hz / 2,500 W	230 V / 50 Hz / 2,500 W	230 V / 50 Hz / 2,500 W
Compressed air	X	–	–
Timer	X	–	–
Water level control	X	X	–
Item no.	10003	10033	10006

# COOLING BATHS

- ✓ Stainless steel bath
- ✓ For 2 or 4 cuvettes
- ✓ Perfect for cooling down the cuvettes



## PRODUCT INFORMATION

	KW 2	KW 4
Dimensions (l x w x h)	350 x 360 x 260 mm	530 x 360 x 260 mm
Water connection	X	X
Item no.	10017	10016



# SR IVOCAP ACCESSORIES

Item no. 20006

## Spring washers for SR Ivocap clamping frame

Regular checking and replacement of the disc springs (approx. every 6 months) guarantees a consistently good work result



Item no. 20075

## Rubber gloves, heat resistant



Item no. 20005

## Floating balls, Ø 20 mm, 1,000 pieces



Item no. 20003

## Dur-Fix-Hardener



Item no. 20146

## Stainless steel lid for IVO-2, IVO-2 special, IVO-2 manual



Item no. 20116

## Stainless steel lid for IVO-4, IVO-4 special, IVO-4 manual



Item no. 20086

## Spiral hose with connector for SR Ivocap-cylinder



# WALL HOLDER


Item no. 20037

Wall holder for 2 SR-Ivocap-cylinder



Item no. 20106

Wall holder for 2 SR-Ivocap-cylinder  
with 3 air sockets




Item no. 20038

Wall holder for 4 SR-Ivocap-cylinder



Item no. 20107


Wall holder for 4 SR-Ivocap-cylinder with 6 air sockets



# STAND HOLDER


Item no. 20073

Stand holder for 2 SR-Ivocap-cylinder




Item no. 20197

Stand holder for 2 SR-Ivocap-cylinder with 3 air sockets




Item no. 20039

Stand holder for 4 SR-Ivocap-cylinder



Item no. 20074

Stand holder for 4 SR-Ivocap-cylinder with 6 air sockets







**BDT GmbH Gesellschaft für Gerätetechnik**  
Industriestraße 27 · D-77656 Offenburg · +49 (0)781 5 86 88 · [info@bdt-dental.de](mailto:info@bdt-dental.de)  
[www.bdt-dental.de](http://www.bdt-dental.de)



as of 01.2022