



PRODUCT GUIDE 2023

Laboratory Edition
January 2023

TERMS

ORDERING

Orders can be placed via authorized Zubler USA personnel during regular business hours; Monday through Friday 8am to 5pm EST, by email to sales@zublerusa.com. Many orders ship the same day if placed by 2pm EST. (except holidays). Correct item numbers are the responsibility of the customer. If you have any questions regarding item numbers or availability, please call customer service for assistance.

TERMS OF SALE

Zubler USA's normal terms of sale are Check, Cash, C.O.D., MasterCard and Visa. Qualifying clients may receive terms of Net 30 upon approval of credit. All shipments invoiced and shipped from Zubler USA must be paid in full within your designated terms. Payments requiring third party fees (such as C.O.D.) will incur additional transaction fees. Credit card transactions over \$500.00 will incur a 3.0% transaction fee.

CREDIT HOLD

Any account with a past due balance may be placed on credit hold and/or incur finance charges at the discretion of Zubler USA until the account and all charges are paid in full.

PRICE CHANGES

Recent cost increases from vendors have given way to multiple price increases. Please be aware that the pricing in this catalog may have changed prior to print. Please note that all prices are subject to change without notice. We advise you to confirm pricing by contacting Zubler USA.

CLAIMS AND RETURNS

"Rejection Period": Buyer is deemed to have accepted the products unless notice of rejection is given within a reasonable period of time, which is agreed to be within (7) days after buyer's receipt of product. Buyer waives any right to invoice rejection or failure of acceptance thereafter. "Return Period": Following the "Rejection Period" there is a 90-day return period from the date of invoice. Any returns of merchandise during either the "Rejection Period" or "Return Period" must be authorized by Zubler USA prior to return to Zubler USA. All returns must be accompanied by a copy of the invoice, a description of the reason for return, as well as any items or information requested by Zubler USA staff. Items returned without the required information will not be considered for credit. All returns are subject to inspection. A 15% restocking fee may be applied to items returned in "Like New" condition. Worn or damaged items will be pro-rated at the discretion of Zubler USA. Credit will not be issued on merchandise returned in open packaging or "poor" condition. No COD Deliveries will be accepted. "Special Order" items or non-stock items such as parts or one time special requests will not be accepted for return or credit

SHIPPING

All shipments will be via UPS Ground unless other arrangements are made at the time the order is placed.

BACK ORDERS

Any item listed as a backorder on an invoice or sales order will ship when the item is available. If you wish to cancel your backorder, please notify customer service.

AGREEMENT TO TERMS AND CONDITIONS


It is understood that you agree to these Terms and Conditions by entering into a business arrangement to purchase products from Zubler USA. We reserve the right to make changes to these Terms and Conditions without notice.

CONTENTS

| | |
|--|------------|
| OVENS | Page |
| Pressing / Ceramic / Sintering | 2 |
| CAD/CAM EQUIPMENT | |
| Milling Machines / Scanners | 6 |
| MILLING MATERIALS | |
| Zirconia / Wax / PMMA / Specialty PMMA | 10 |
| CERAMICS | |
| 9.2 Ceramic / ConceptArt Stains / ConceptPress Lithium Disilicate | 13 |
| INVESTING & PRESSING | |
| Materials / Equipment / Accessories | 16 |
| HAND PIECES & ROTARY | |
| Hand Pieces / Rotary Polishers and Burs | 18 |
| SUCTION | |
| Central Suction / CadCam Suction/ Workbench Suction/ Accessories | 22 |
| REMOVABLE EQUIPMENT & LAB ESSENTIALS | |
| Boil Out / Polymerization / Duplicating Material / Abrasives / Putty | 27 |
| TERMS | LEFT |
| CONTACT INFORMATION & LOCATIONS | BACK COVER |


ZUBLER VARIO OVENS

Sintering, Ceramic, Pressing & Accessories




MADE IN GERMANY

| OVEN MODEL | VARIO S430 | VARIO230ZR | VARIOPRESS300.Ezr |
|------------|------------------------------------|--------------------------------------|---|
| TYPE | Zirconia Sintering | Ceramic & Glazing | Pressing, Ceramic & Glazing |
| FEATURES | TTC Linear Colling/Speed Sintering | TTC Linear Cooling / Z120 Technology | Fully Automatic Pressing / TTC Linear Cooling |
| ITEM NO. | 898/430 | 898/046 | 898/038 |
| PRICE | \$12,950.00 | \$6,850.00 | \$9,190.00 |
| WARRANTY | 1 YR | 2 YR | 2 YR |
| NOTES: | Includes 2 trays & beads | Price Includes Vacuum Pump | Price includes Vacuum Pump |



Zubler's quiet yet powerful P3 vacuum pump is included with the purchase of Vario ceramic and pressing ovens. Designed for quiet operation, yet powerful vacuum.

| | |
|----------|----------|
| Item No: | Price: |
| 898/053 | \$850.00 |



ZUBLER VARIO - OVEN PARTS & ACCESSORIES

| PIC# | ITEM DESCRIPTION | ITEM# | PRICE |
|------|---|----------|----------|
| 1 | Honeycomb Firing Tray with Pegs | 898/1045 | \$74.99 |
| 2 | Ceramic Firing Insulation (V200/V230) | 898/105 | \$69.99 |
| 3 | Ceramic Firing Insulation (V300/VP300/VP330) | 898/110 | \$69.99 |
| 4 | Pressing Insulation Tray Insert (VP300/VP330) | 898/109 | \$49.99 |
| 5 | Pressing Insulation Tray (VP300/VP330) | 898/108 | \$69.99 |
| 6 | Pressing Insulation Tray W/Insert (VP300/VP330) | 898/107 | \$114.99 |
| - | Fuse F1 (15 Amp for V200/V230/V300/VP300/VP330) | 028/0328 | \$3.85 |
| - | Fuse F3 (5 Amp for V200/V230/V300/VP300/VP330) | 028/0324 | \$3.85 |

| PIC# | ITEM DESCRIPTION | ITEM# | PRICE |
|------|---|----------|----------|
| 7 | Zubler Tongs | 898/4136 | \$162.99 |
| 8 | Zubler Tweezers | 898/106 | \$35.00 |
| 9 | Sintering Dish - Lid only (S400/S430) | 896/0421 | \$142.00 |
| 9 | Sintering Dish - Dish only (S400/S430) | 896/0420 | \$204.00 |
| 10 | Sintering Dish Support Ring (S400/S430) | 896/0406 | \$92.00 |
| 11 | Sintering Tray Insulation Level 1 (S400/S430) | 896/0408 | \$141.99 |
| 12 | Sintering Tray Insulation Level 2 (S400/S430) | 896/0407 | \$149.00 |
| - | Fuse (16 Amp for S400/S430) | 028/0129 | \$5.35 |

MIHM VOGT
Sintering Ovens & Accessories



MV-R


TABEO-1/M/ZIRKON-100

HTS-2/M/ZIRKON-120

HT-2/M/ZIRKON-120

| | | | |
|---|--|---------------------------------------|--------------------------------------|
| Zirconia Sintering | Zirconia Sintering | Zirconia Sintering | Zirconia Sintering |
| 40 Single Crowns- | 20 Single Crowns | 60 Single Crowns | 90 Single Crowns |
| Two 10cm tray capacity | Single 10cm Tray Capacity | Two 12cm tray capacity | Three 12cm Tray Capacity |
| Touch Screen / Speed Sintering / Rapid Sintering up to 200C/Min Heat Rate/ USB Port | 25°C/Min Heat Rate / Regular Sintering | 120°C/Min Heat Rate / Speed Sintering | 99°C/Min Heat Rate / Speed Sintering |
| 7204000100 | 00009270 | 00009550 | 00009650 |
| \$ 14,700.00 | \$ 8,890.00 | \$ 13,500.00 | \$ 15,400.00 |





Mihm Vogt Sintering Oven Parts & Accessories

NOTE: * = For items to be used with classic sintering cycles with heat/cool rates less than 25°C/min.

| PIC# | ITEM DESCRIPTION | SIZE | ITEM# | PRICE |
|------|---------------------------------|-------|------------|--------------|
| 1 | Sintering Dish Lid (*) | 120mm | 00096025 | \$199.00 |
| | | 100mm | 00095025 | \$176.00 |
| 2 | SiC Sintering Tray (NEW!) | 120mm | 8201000047 | \$685.00 |
| | | 100mm | 8201000048 | \$613.00 |
| 3 | MV Thermal Fabric (Pkg of 5) | 120mm | 8203000033 | \$149.00 |
| 4 | MV Thermal Fabric (Pkg of 5) | 40mm | 8203000039 | \$85.00 |
| 5 | Sintering Dish (*) Unventilated | 120mm | 07100002 | Discontinued |
| | | 100mm | 07100001 | Discontinued |
| 6 | Speed Sintering Dish (Slotted) | 120mm | 07100006 | \$213.00 |
| 7 | Mihm Vogt Fan Battery | | 00096035 | \$107.00 |
| 8 | M/V Power Cord | | 00096042 | \$40.00 |
| 9 | Ring of Speed Sintering Tray | 120mm | 00096038 | \$232.00 |
| | | 100mm | 00095038 | \$215.00 |
| 10 | Speed Sintering Tray | 120mm | 00096037 | \$249.00 |
| | | 100mm | 00095037 | \$231.00 |

| PIC# | ITEM DESCRIPTION | SIZE/FITMENT | ITEM# | PRICE |
|------|-------------------------------|---------------------------------|------------|----------|
| 11 | Sinter Granulate | 200g | 00095023 | \$153.00 |
| 12 | Sintering Beads (2mm) | 200g | 00096023 | \$130.99 |
| 13 | Tray Support Disc | HT-2M/HTS-2M | 7300006 | \$177.00 |
| 14 | Upper Door Insulation | HT-2M/HTS-2M | 7300003 | \$193.00 |
| | | HT (all)/HTS-2M | 00096078 | \$195.99 |
| 15 | Lower Door Insulation | HT-S/HTS-SPEED | 00095078 | \$239.99 |
| | | MV-R (Not Pictured) | 8203000038 | \$178.99 |
| 16 | MV Heating Elements | Call us for pricing and fitment | | \$Call |
| 17 | Tray Stability Pins | HT / HTS / SPEED | 00096080 | \$33.99 |
| N/A | Door Insulation Stability Pin | MV-R | 8203000037 | \$15.99 |
| 18 | Upper Door Insulation | HT / HT-SPEED, | 00096079 | \$169.00 |
| | | HTS / HTS-SPEED | 00095079 | \$183.00 |
| 19 | Bowl Support Pins (W/Flair) | MV-R (Not Pictured) | 8203000041 | \$206.99 |
| | | HT / HTS / SPEED | 00096081 | \$124.00 |
| 20 | Honeycomb Cooling Tray | All | 995415 | \$65.00 |

Important!: When using SiC Dishes, you must place MV Thermal Fabric under and between each dish. Thermal Fabric may also be used as a lid.

DCS MILLING MACHINES

5 Axis, Simultaneous, Wet/Dry Milling



DCS - DC 1

| | |
|-------------------|---|
| DESCRIPTION: | 5 Axis Simultaneous Wet/Dry mill |
| Pricing: | Contact Zubler USA for Quotation |
| Milling Range: | 360° tilt angle (B-axis) ± 30° milling angle (A- / B-axis) |
| Max Speed: | 60,000 rpm |
| Spindle: | High Frequency, Ball Bearing |
| Tool Size: | 6mm |
| Bur Capacity: | 18 |
| Bur Changing: | Automatic |
| Disc Capacity: | 1 |
| Accuracy: | ± 0.01 mm (10 Microns) |
| Calibration: | Automatic |
| Dimensions D/W/H: | 87.5mm x 55mm x 93.5mm |
| Materials: | Zirconia, Wax, PMMA, Titanium, CoCr, Wood, Resin, Acrylic |
| Compressed Air: | 7 Bar |
| Voltage: | 230Vac |

DCS - DC 7

| | |
|-------------------|---|
| DESCRIPTION: | 5 Axis Simultaneous Wet/Dry mill |
| Pricing: | Contact Zubler USA for Quotation |
| Milling Range: | 360° tilt angle (B-axis) ± 30° milling angle (A- / B-axis) |
| Max Speed: | 60,000 rpm |
| Spindle: | High Frequency, Ball Bearing |
| Tool Size: | 6mm |
| Bur Capacity: | 20 (Exchangable Mag) |
| Bur Changing: | Automatic |
| Disc Capacity: | 7 |
| Accuracy: | ± 0.01 mm (10 Microns) |
| Calibration: | Automatic |
| Dimensions D/W/H: | 73cm x 195cm x 81.5cm |
| Materials: | Zirconia, Wax, PMMA, Titanium, CoCr, Wood, Resin, Acrylic |
| Compressed Air: | 7 Bar |
| Voltage: | 230Vac |

ZUBLER MILLING MACHINES

New Technology and Robust Capabilities



Z- 565 DC

Z-515

Z-505

DESCRIPTION:

Pricing:

Milling Range:

Max Speed:

Spindle:

Tool Size:

Bur Capacity:

Bur Changing:

Disc Capacity:

Accuracy:

Calibration:

Dimensions D/W/H:

Materials:

Compressed Air:

Voltage:

12 Disc C-Clamp Auto Changing 5 Axis Simultaneous Dry mill

Contact Zubler USA for Quotation

X:120mm/Y:150mm/Z:80mm /
A: ±35° / B: 360°

60,000rpm

1.5kW Water Cooled

4mm

8

Automatic

1

± 0.01 mm (10 Microns)

Automatic

820mm x 750mm x 720mm

Zirconia, Wax, PMMA

6 Bar

220Vac

5 Axis Simultaneous Dry mill with C-Clamp

Contact Zubler USA for Quotation

X:120mm/Y:150mm/Z:80mm /
A: 360°

60,000rpm

1.5kW Water Cooled

4mm

3

Automatic

1

± 0.01 mm (10 Microns)

Automatic

820mm x 750mm x 720mm

Zirconia, Wax, PMMA

6 Bar

220Vac

4 Axis Wet mill

Contact Zubler USA for Quotation

X:125mm/Y:130mm/Z:80mm /
A: 360°

60,000rpm

1.5kW Water Cooled

4mm

6

Automatic

1

± 0.01 mm (10 Microns)

Automatic

730mm x 720mm x 610mm

Glass Ceramic, PMMA, Resin

6 Bar

220Vac

MEDIT T-SERIES SCANNERS

Where accuracy begins...



Contact us for a quotation

Super Fast Scanning!
Extreme reality!



| MEDIT Desktop Scanners | T710 | T510 | T310 |
|------------------------|--------------------------------------|----------------------|----------------------|
| Resolution of Camera | Mono 5.0 (MP) x 4 | Mono 5.0 (MP) x 2 | Mono 5.0 (MP) x 2 |
| Point spacing | 0.040 mm | | |
| Scan Area | 100mm x 73mm x 60mm | | |
| Scan Principle | Phase-shifting optical triangulation | | |
| Size | 505mm x 271mm x 340 mm | | |
| Weight | 15 kg | | |
| Light Source | LED, 150 ANSI-lumens | | |
| Light source | Blue LED | | |
| Connection | USB 3.0 B Type | | |
| Power | AC 100-240V, 50-60 Hz | | |
| Accuracy (ISO 12836) | 4µm | 7µm | 9µm |
| Full arch scan speed | 8 sec(8cut) | 12 sec(8cut) | 18 sec(8cut) |

| MEDIT Scanners & Modules | T710 | T510 | T310 |
|--------------------------|----------|----------|----------|
| Base Price | \$12,900 | \$9,000 | \$6,000 |
| Color texture | included | included | included |
| Impression scanning | included | X | X |
| Flexible scanning | included | included | \$1,100 |
| Articulator scanning | included | included | \$1,100 |
| Replica denture | included | \$1,100 | \$1,100 |
| Orthodontic scanning | included | \$1,100 | \$1,100 |
| Total | \$12,900 | \$11,200 | \$10,400 |




Unsurpassed speed, accuracy and reliability - Rediscover your productivity.



INTRAORAL SCANNERS

Medit i600, i700 & i700 Wireless Scanners



Fast, accurate scans, from a lightweight and comfortable scanner packed with features!
Scanning Reimagined!

MEDIT Intraoral Scanners

| | |
|---------------------|---|
| Scanning frame | Up to 70 FPS |
| Imaging technology | 3D-in-motion video 3D full color capture |
| Light source | LED |
| Anti-fogging | Adaptive anti-fogging |
| Full-arch | 10.9µm ± 0.98 |
| Dimensions | 248 x 44 x 47.4 mm |
| Weight | 245 g |
| Tip size | 22.2 x 15.9 mm |
| Mirror angle | 45-degree angle |
| Scan area | 15 x 13 mm |
| Reversible tip | Yes |
| Remote control mode | Yes |
| UV-C disinfection | Yes |


MEDIT i700

| | |
|--|------------|
| Calibration Tool | \$350.00 |
| Tip set(4ea) | \$300.00 |
| Desktop Cradle | \$60.00 |
| Wall Mount Holder | \$35.00 |
| Medical Adapter & Power Cord | \$50.00 |
| Power Hub | \$130.00 |
| Power Delivery Cable (2.0m) | \$150.00 |
| USB 3.0 Cable | \$50.00 |
| Wrist Strap | \$15.00 |
| i700 Handpiece Cover | \$10.00 |
| Extended warranty (2 additional years) | \$1,700.00 |

\$Contact Us

MEDIT i700 Wireless

\$Contact Us.



ZIRCONIA DISCS

Sinter full arches without support structures faster than you ever believed was possible!



ZUBLER 3:45 Full Arch Speed Zirconia

| ITEM# | SHADE | SIZE (mm) | PRICE |
|------------|-------|-------------|----------|
| W1100A025 | WHITE | 98mm x 25mm | \$145.00 |
| W1100A030 | WHITE | 98mm x 30mm | \$160.00 |
| W1100A035 | WHITE | 98mm x 35mm | \$173.00 |
| W1100AB125 | B1 | 98mm x 25mm | \$145.00 |
| W1100AB130 | B1 | 98mm x 30mm | \$160.00 |
| W1100BL125 | BL1 | 98mm x 25mm | \$145.00 |
| W1100BL130 | BL1 | 98mm x 30mm | \$160.00 |
| W1100BL325 | BL3 | 98mm x 25mm | \$145.00 |
| W1100BL330 | BL3 | 98mm x 30mm | \$160.00 |
| W1100XA125 | A1 | 98mm x 25mm | \$145.00 |
| W1100XA130 | A1 | 98mm x 30mm | \$160.00 |

More Shades and sizes coming soon and may already be available, get in touch with us to find out!

Color liquids also available.

ZUBLER 3:45 Technical Specifications

| MATERIAL | UNIT | 3:45 ZIRCONIA | METHOD |
|-------------------------------|-------|---------------|----------------------------|
| ZrO2 | Wt. % | 88~90 | ICP |
| Y2O3 | | 7.0-8.0 | |
| SiO2 | | ≤ 0.01 | |
| Fe2O3 | | ≤ 0.001 | |
| CaO | | ≤ 0.007 | |
| Na2O | | ≤ 0.004 | |
| Property | Unit | 345 Zirconia | Test Method |
| Sintered Density | g/cm3 | ≥6.04 | Archimedes' Method |
| Bending Strength | MPa | 1100±50 | 3-Point-Bending ISO 6872 |
| Transparency | % | ≥ 46 | Spectrophotometer (T=p/p0) |
| Thermal Expansion Coefficient | K-1 | 10.5 - 10.6 | ASTM Test Method E 289 |

MILLING WAX Wax

| | | | |
|------|----------|-------------|---------|
| Grey | 898/9001 | 98mm x 20mm | \$35.00 |
|------|----------|-------------|---------|

The perfect milling wax, beautiful marginal integrity, meets easy handling.

PMMA
Mono Shaded and Denture



MONO SHADED PMMA

| SHADE | 25mm | PRICE | 30mm | PRICE | 35mm | PRICE |
|-------|----------|---------|----------|---------|---------|---------|
| A1 | 2525A1 | \$36.00 | 2530A1 | \$45.00 | 2535A1 | \$54.00 |
| A2 | 2525A2 | \$36.00 | 2530A2 | \$45.00 | 2535A2 | \$54.00 |
| A3 | 2525A3 | \$36.00 | 2530A3 | \$45.00 | | |
| A3.5 | 2525A3.5 | \$36.00 | 2530A3.5 | \$45.00 | | |
| A4 | 2525A4 | \$36.00 | 2530A4 | \$45.00 | | |
| B1 | 2525B1 | \$36.00 | 2530B1 | \$45.00 | 2535B1 | \$54.00 |
| B2 | 2525B2 | \$36.00 | 2530B2 | \$45.00 | | |
| B3 | 2525B3 | \$36.00 | 2530B3 | \$45.00 | | |
| B4 | 2525B4 | \$36.00 | 2530B4 | \$45.00 | | |
| BL2 | 2525BL2 | \$36.00 | 2530BL2 | \$45.00 | | |
| BL3 | 2525BL3 | \$36.00 | 2530BL3 | \$45.00 | 2535BL3 | \$54.00 |
| C1 | 2525C1 | \$36.00 | 2530C1 | \$45.00 | | |
| C2 | 2525C2 | \$36.00 | 2530C2 | \$45.00 | | |

Note: Multilayered discs are available also!
See the Zubler USA - Newstetic Catalog for Details.

MONO SHADED PMMA (Continued)

| SHADE | 25mm | PRICE | 30mm | PRICE |
|-------|--------|---------|--------|---------|
| C3 | 2525C3 | \$36.00 | 2530C3 | \$45.00 |
| C4 | 2525C4 | \$36.00 | | |
| D2 | 2525D2 | \$36.00 | 2530D2 | \$45.00 |
| D3 | 2525D3 | \$36.00 | 2530D3 | \$45.00 |
| D4 | 2525D4 | \$36.00 | 2530D4 | \$45.00 |

DENTURE PMMA

| SHADE | 25mm | PRICE | 30mm | PRICE |
|-------------|----------------|---------|----------------|---------|
| PINK-Veined | 2525PINK | \$36.00 | 2530PINK | \$40.00 |
| A1/PINK | 25PINK-ML25A1 | \$85.00 | 25PINK-ML30A1 | \$90.00 |
| A2/PINK | 25PINK-ML25A2 | \$85.00 | 25PINK-ML30A2 | \$90.00 |
| B1/PINK | 25PINK-ML25B1 | \$85.00 | 25PINK-ML30B1 | \$90.00 |
| BL2/PINK | 25PINK-ML25BL2 | \$85.00 | 25PINK-ML30BL2 | \$90.00 |
| BL3/PINK | 25PINK-ML25BL3 | \$85.00 | 25PINK-ML30BL3 | \$90.00 |

PMMA MULTILAYER & CAD/CAM

Process Accessories




PORTUX MULTILAYER PMMA


| PRICE -> | \$76.21 | \$77.09 | \$80.70 | \$87.93 | \$89.13 |
|----------|----------|----------|----------|----------|----------|
| SHADE | 18mm | 20mm | 25mm | 30mm | 35mm |
| A1 | NW 21846 | NW 21858 | NW 22145 | NW 23116 | NW 22576 |
| A2 | NW 21849 | NW 21861 | NW 22148 | NW 23117 | NW 24457 |
| A3 | NW 22083 | NW 22078 | NW 22151 | NW 23118 | |
| A3.5 | | | NW 22225 | NW 24671 | |
| B1 | NW 21852 | NW 21864 | NW 22154 | NW 23119 | NW 22577 |
| B2 | NW 22350 | NW 22355 | NW 22360 | NW 23120 | |
| B3 | | | NW 26921 | NW 27157 | |
| BL2 | NW 21954 | NW 21933 | NW 22332 | NW 24141 | NW 24147 |
| BL3 | NW 27144 | NW 27145 | NW 27146 | NW 27147 | NW 27148 |
| C1 | NW 24476 | NW 24477 | NW 24478 | NW 24868 | |
| C2 | NW 24479 | NW 24480 | NW 24481 | NW 24870 | |
| C3 | | | NW 26922 | NW 27160 | |
| D2 | NW 24482 | NW 24483 | NW 24484 | NW 24872 | |
| D3 | NW 23336 | NW 23337 | NW 23338 | NW 23339 | |

CAD/CAM PROCESS ACCESSORIES

| DESCRIPTION | ITEM NO. | PRICE |
|--------------------------------|----------|---------|
| DC Resin Elite Translucent 8gr | 995115 | \$95.99 |
| DC Resin Elite Opaque 8gr | 995116 | \$95.99 |
| DC Resin Elite White 8gr | 995117 | \$95.99 |
| Clear Scan Powder W/Brush | 995425 | 49.95 |



DC CERAM CONCEPT PRESS Lithium Disilicate



Lifelike Fluorescence! Stackable 2g & 3g Ingot sizes!

DENTIN INGOTS (Medium Translucency)

| SHADE | ITEM# | SIZE | PRICE | ITEM# | SIZE | PRICE |
|-------|---------|------|---------|---------|------|----------|
| A1 | 47-0212 | 2g | \$95.00 | 47-0012 | 3g | \$115.00 |
| A2 | 47-0213 | 2g | \$95.00 | 47-0013 | 3g | \$115.00 |
| A3 | 47-0214 | 2g | \$95.00 | 47-0014 | 3g | \$115.00 |
| A3.5 | 47-0215 | 2g | \$95.00 | 47-0015 | 3g | \$115.00 |
| A4 | 47-0216 | 2g | \$95.00 | 47-0016 | 3g | \$115.00 |
| B1 | 47-0222 | 2g | \$95.00 | 47-0022 | 3g | \$115.00 |
| B2 | 47-0223 | 2g | \$95.00 | 47-0023 | 3g | \$115.00 |
| B3 | 47-0224 | 2g | \$95.00 | 47-0024 | 3g | \$115.00 |
| B4 | 47-0225 | 2g | \$95.00 | 47-0025 | 3g | \$115.00 |



COLOR TRANSPARENT Ingots (Hi Translucency)

| SHADE | ITEM# | PRICE | ITEM# | PRICE |
|-------|---------|------------|---------|-------------|
| CT1 | 47-4211 | 2g \$95.00 | 47-4011 | 3g \$115.00 |
| CT2 | 47-4212 | 2g \$95.00 | 47-4012 | 3g \$115.00 |
| CT3 | 47-4213 | 2g \$95.00 | 47-4013 | 3g \$115.00 |

INTENSIVE DENTIN Ingots (Hi Opacity)


| ITEM# | PRICE | ITEM# | PRICE |
|-------|--------------------|---------|-------------|
| ID1 | 47-3212 2g \$95.00 | 47-3012 | 3g \$115.00 |
| ID2 | 47-3213 2g \$95.00 | 47-3013 | 3g \$115.00 |
| ID3 | 47-3214 2g \$95.00 | 47-3014 | 3g \$115.00 |
| ID4 | 47-3215 2g \$95.00 | 47-3015 | 3g \$115.00 |
| ID5 | 47-3216 2g \$95.00 | 47-3016 | 3g \$115.00 |

PEARL Ingots (Opalescent)

| ITEM# | PRICE | ITEM# | PRICE |
|-------|--------------------|---------|-------------|
| PRL1 | 47-4221 2g \$95.00 | 47-4021 | 3g \$115.00 |
| PRL2 | 47-4222 2g \$95.00 | 47-4022 | 3g \$115.00 |
| PRL3 | 47-4223 2g \$95.00 | 47-4023 | 3g \$115.00 |



DC CERAM 9.2
Ceramic System



DC CERAM 9.2 Ceramic Powders

| Description: | Available in: | Size: | Price: |
|-------------------------------|--|-------|----------|
| Liner Powders | 1, 2, 3, 4, 5, N | 15g | \$40.99 |
| Chroma Dentin Powders | Shades A1 - D4, SLA0, SLB0 | 20g | \$48.99 |
| Dentin Modifiers | Shades A, B, C, D, White, Yellow, Orange, Brown, Pink, Violet, Blue | 20g | \$47.99 |
| Dentin Powders | Shades A1 - D4, SLA0, SLB0, Dentin Correction | 20g | \$47.99 |
| Special Dentin 3 - 16 | Bright, Medium, Dark | 20g | \$47.99 |
| Special Dentin 3 - 16 (Large) | Bright, Medium, Dark | 70g | \$142.00 |
| Incisal Powders | Enamel 1, 2, 3, 4, Neutral White, Correction Material | 20g | \$47.99 |
| Gingival (Pink) Powders | Gingiva 1, Gingiva 2 | 20g | \$47.99 |
| Margin Powders | 1, 2, 3, 4, 5, 6 | 20g | \$57.99 |
| Opal Powders | Opal 1, Opal 2 | 20g | \$61.99 |
| Transparent Powders | Neutral, Clear, Yellow, Orange, Brown, Pink, Light Blue, Blue, Grey, White | 20g | \$55.99 |

DC CERAM 9.2 Liquids & Shade Guides


| Liquid | Size: | Price: |
|-----------------|-------|---------|
| Liner Liquid | 100ml | \$36.99 |
| Build Up Liquid | 100ml | \$36.99 |
| Stain / Glaze | 25ml | \$21.99 |
| Stain / Glaze | 100ml | \$71.00 |

| 9.2 Shade Guides | Price: |
|----------------------------|----------|
| Enamel/Opal/Transpa | \$140.00 |
| Modifier / Gingiva / Liner | \$140.00 |

Note: Individual item numbers for ceramics and related components are too numerous to print. Provide us with a description of the system, component and shade or application and size along with the quantity when placing your order.



DC CERAM CONCEPT ART
Fluorescent Stains



CONCEPT ART STAINS Flourcent Shades and Colors

| Description | Item No. | Size | Price |
|------------------------|----------|------|---------|
| Gum 4 Flame Red | 41-1035 | 3g | \$59.50 |
| Black Fluorescent | 41-1024 | 3g | \$59.50 |
| Grey Blue Fluorescent | 41-1027 | 3g | \$59.50 |
| Dark Brown Fluorescent | 41-1026 | 3g | \$59.50 |
| Yellow Fluorescent | 41-1016 | 3g | \$59.50 |
| White Fluorescent | 41-1015 | 3g | \$59.50 |
| Moss Fluorescent | 41-1022 | 3g | \$59.50 |
| Violet Fluorescent | 41-1019 | 3g | \$59.50 |

| Description | Item No. | Size | Price |
|---------------------|----------|------|---------|
| Shade C Fluorescent | 41-1012 | 3g | \$59.50 |
| Orange Fluorescent | 41-1017 | 3g | \$59.50 |
| Shade D Fluorescent | 41-1013 | 3g | \$59.50 |
| Khaki Fluorescent | 41-1021 | 3g | \$59.50 |
| Smoke Fluorescent | 41-1025 | 3g | \$59.50 |
| Shade B Fluorescent | 41-1011 | 3g | \$59.50 |
| Shade A Fluorescent | 41-1010 | 3g | \$59.50 |
| Vanilla Fluorescent | 41-1028 | 3g | \$59.50 |



Fluorescent stains

Compared to classic stains, fluorescent stains have a huge advantage. The fluorescent components operate to amplify light. This means the concept Art stains will retain their light optical appearance (similar to that of a natural tooth) even in relatively dark areas of the mouth cavity. Another advantage to the concept Art fluorescent stains is that graying does not occur.

Basic Kit



Professional Kit



FLEX RING FORMERS
Longer lasting with extrusions



| FLEX Ring (Red Ring Formers) | Item No. | Price |
|------------------------------|----------|---------|
| For 100 gram rings | 558/150 | \$33.00 |
| For 200 gram rings | 558/1512 | \$36.00 |
| For 300 gram rings | 558/1513 | \$58.00 |


| FLEX Ring Level Tops | Item No. | Price |
|----------------------|------------|---------|
| For 100 gram rings | 563/075101 | \$8.50 |
| For 200 gram rings | 563/07502 | \$10.25 |
| For 300 gram rings | 563/07505 | \$11.75 |

| FLEX Ring Bases | Item No. | Price |
|-----------------------------------|-----------|---------|
| 100 gram - Ingot Diameter (12 mm) | 563/07508 | \$15.00 |
| 200 gram - Ingot Diameter (12 mm) | 563/07501 | \$15.00 |
| 200 gram - Ingot Diameter (16 mm) | 563/07530 | \$35.00 |
| 300 gram - Ingot Diameter (12 mm) | 563/07503 | \$35.00 |
| 300 gram - Ingot Diameter (16 mm) | 563/07504 | \$39.00 |

| FLEX Ring Base and Top - 3 Pack | Item No. | Price |
|---------------------------------|----------|---------|
| For 100 gram rings (Grey 13mm) | 898/118 | \$31.00 |
| For 200 gram rings (Grey 13mm) | 898/119 | \$33.00 |



INVEST & PRESS Materials and Equipment



Zubler Investment

| Material | Size / Qty. | Item No. | Price |
|---------------------------------|-------------|----------|----------|
| ZUBLER 144 Universal Speed Inv. | 100 x 100g | 898/2035 | \$124.95 |
| CONCEPTVEST LiSi2 Pressing Inv. | 50 x 100g | 898/2065 | \$85.00 |
| ZUBLER 144 Investment Liquid | 1 Liter | 898/204 | \$18.95 |

ZUBLER Disposable Plungers

| | Item No. | Price |
|---------------------------------------|----------|---------|
| Ø 12mm (Leucite Ceramics) | 20/0235 | \$56.95 |
| Ø 13mm (For lithium disilicate) | 20/0238 | \$56.95 |
| Ø 13mm - For e.max® (*) "Press-Multi" | 20/0241 | \$56.95 |
| Ø 16mm - For large & 300g rings | 20/0237 | \$47.00 |

Accessories for Investing, Pressing & Casting

| | Item No. | Price |
|-------------------------------|----------|----------|
| ZUBLER tongs | 898/4136 | \$162.99 |
| ZUBLER Tweezers (Length 30cm) | 898/106 | \$35.00 |

VARIO Balance Dosing Device

898/420

\$425.00

Dosing for Investment, Gypsum and Acrylic. 100% predictable results every time!

VARIO Surface Conditioner

898/425

\$509.99

Improves surfaces of pressed ceramics, prevents craters and other surface anomalies. Used with hot rings immediately from the press oven.



HIGH SPEED HAND PIECES

From Apple Dental



APPLE Dental High Speed Hanpieces

| Item No. | Description | Price |
|-----------------|---|----------|
| A1-TU-M4 | Apple A1 Handpiece | \$110.00 |
| A2-QC-M4 | Apple A2 HP w/ LED & Quick Connect | \$170.00 |
| A2-TUP-M4 | Apple A2 Handpiece with LED | \$140.00 |
| Apple A45-M4 | Apple 45° Surgical Handpiece | \$185.00 |
| Apple Gold-L-M4 | Apple Gold-L-M4 Slow Speed Kit (Straight / Contra Angle / Quick Connect) | \$250.00 |
| Apple-SH HP | Apple-SH Slow Speed HP | \$110.00 |
| Apple-AM-M4 | Apple AM-M4 Motor Only | \$75.00 |
| A180205 | Apple Portable FC Air Regulator w/FT Pdl | \$165.00 |



VARIO STAR HANDPIECES

From Zubler GmbH



VARIO Star T50

Bench-top System

| | |
|---------------------|---|
| Included Equipment: | Control Unit, Handpiece, Cradle, Foot Pedal |
| Optional Equipment: | Knee Pedal |
| Item No. | 804/3151 |
| Price: | \$1,195.00 |

Accessories:

| Item No. | Price |
|----------|----------|
| 896/1021 | \$119.00 |
| 896/1016 | \$99.00 |
| 896/1018 | \$15.00 |
| 896/1017 | \$14.00 |
| 804/1017 | \$205.00 |
| 804/1016 | \$210.00 |
| 896/1021 | \$119.00 |



VARIO Star K50

Knee Control System

| | |
|---------------------|--------------------------------------|
| Included Equipment: | Knee Control Unit, Handpiece, Cradle |
| Optional Equipment: | Active Cradle |
| Item No. | 804/3111 |
| Price: | \$1,245.00 |

Polishers and wheels on the next page!

HANDPIECE POLISHERS AND ABRASIVES

For low speed handpiece systems

- AWG/BWG – (Shown left) Adjusting Un-sintered Zirconia
- DGO - Adjustment of Zirconia & Ceramic
- TCP – PANTHER Zirconia, Ceramics & Abutments
- PSC – PANTHER Diamond Grinding Stones
- EZA – Finishing Sintered Zirconia
- CCE – Ceramic Polisher
- MCE – Ceramic Adjustment
- SCE – Ceramic Polishers



GREEN STATE TOOLS

| Item No. | Shape | Color | Grit | Material Usage | Ø (mm) | [H] (mm) | Rpm Opt/Max | Price |
|-------------------|------------|-------|----------|-------------------------------------|--------|----------|-------------|-------------|
| AWGBWGSET | SET | | | Med/Fine Unsintered Zirconia Kit | | | | \$18.50 Ea. |
| BWG164104MF | XL Flame | Green | Med/F | Medium Fine for unsintered Zirconia | 10 | 22.7 | 8k/20k | \$6.25 Ea. |
| AWG106104F | Flame | White | Fine | Unsintered Zirconia | 4.3 | 12 | 8k/20k | \$1.95 Ea. |
| AWG111104M | Flame | Green | Med | Unsintered Zirconia | 5.6 | 13 | 6k/20k | \$2.05 Ea. |
| AWG115104M | LG Flame | Green | Med | Unsintered Zirconia | 7 | 13 | 6k/20k | \$2.05 Ea. |
| AWG146104M | Knife Edge | Green | Med | Unsintered Zirconia | 15 | 0.3 | 6k/20k | \$2.20 Ea. |
| AWG176104M | Cylinder | Green | Med | Unsintered Zirconia | 6 | 15 | 6k | \$3.50 Ea. |
| AWG560104MF | Cylinder | White | Med/F | Unsintered Zirconia | 6 | 15 | 6k | \$5.20 Ea. |
| Jay's GS Tool Kit | Kit | | Multiple | Unsintered "Green State" Zirconia | | | | \$89.00 Ea. |


Jay's GS Tool Kit - Eliminate post sintering processing! More information and Jay's instructional videos can be found here: <https://zublerusa.com/gs-tool-kit/>



Jay Emirzian CDT

Jay's GS Tool Kit may be purchased as a complete set or individually. Individual item information is listed here:

<https://zublerusa.com/gs-tool-kit/>



AWGBWGSET



BWG164104MF




AWG560104MF



AWG115104M



AWG146104M



AWG176104M

HANDPIECE POLISHERS AND ABRASIVES

For low speed handpiece systems




Green Diamond Grinders


| Item No. | Shape | Color | Material Usage | Ø (mm) | [H] (mm) | Price |
|----------|-----------|-------|--------------------|--------|----------|------------|
| DG001 | Cylinder | Green | Zirconia & Ceramic | 12 | 0.5 | \$12.00 Ea |
| DG002 | Trapazoid | Green | Zirconia & Ceramic | 13 | 0.5 | \$12.00 Ea |
| DG012 | Wheel SM | Green | Zirconia & Ceramic | 2 | 1.2 | \$18.00 Ea |
| DG007 | Wheel MD | Green | Zirconia & Ceramic | 3.5 | 1.5 | \$15.00 Ea |
| DG004A | Lense SM | Green | Zirconia & Ceramic | 1.6 | 0.7 | \$12.00 Ea |
| DG004 | Lense LG | Green | Zirconia & Ceramic | 2.2 | 0.7 | \$15.00 Ea |
| DG005 | Wheel LG | Green | Zirconia & Ceramic | 22 | 2 | \$16.00 Ea |



DG001



DG002



DG012



DG007



DG004A



DG004





DG005

Premium green diamond stone grinders for sintered zirconia and all ceramics. Exceptional heat reduction & tool life!

HANDPIECE POLISHERS AND ABRASIVES


For low speed handpiece systems



EZA105104XC




EZA140104XC




EZA177104XC



EZA200104XC




SCE100104CGE1Z



SCE235104CGE1R



MCE138104C




MCE191104



MCE210104C



MCE239104C



MCE239900M



CCE235900M

TOOLS FOR SINTEREED ZIRCONIA, CERAMIC AND METAL

| Item No. | Shape | Color | Grit | Material Usage | Ø (mm) | [H] (mm) | Rpm Opt/Max | Price |
|----------------|------------------|---------|---------|--|--------|----------|-------------|-------------|
| EZA105104XC | Flame | Green | XtraCrs | Sintered Zirconia (Suitable for Ceramics & TITN) | 5.5 | 16 | 15k/20k | \$8.25 Ea. |
| EZA140104XC | Disc | Green | XtraCrs | Sintered Zirconia (Suitable for Ceramics & TITN) | 26.2 | 2 | 6k/20k | \$40.50 Ea. |
| EZA177104XC | Cylinder | Green | XtraCrs | Sintered Zirconia (Suitable for Ceramics & TITN) | 6 | 15 | 15k/20k | \$8.75 Ea. |
| EZA200104XC | Whl Flat Edge | Green | XtraCrs | Sintered Zirconia (Suitable for Ceramics & TITN) | 22 | 3.2 | 6k/20k | \$20.50 Ea. |
| SCE100104CGE1Z | Flame | Blue | Fine | Ceramic (Suitable for ZIRC, NP & TITN) | 5.5 | 18 | 6k/20k | \$9.95 3Pk |
| SCE235104CGE1R | Lense Knife edge | Blue | Fine | Ceramic (Suitable for ZIRC, NP & TITN) | .22 | 3.5 | 6k/20k | \$22.99 Ea |
| MCE138104C | Cylinder | Yellow | Crse | Sintered Zirconia & Ceramic (Suitable for NP & TITN) | 4 | 10.5 | 15k/20k | \$17.95 Ea. |
| MCE191104M | Cup | Seafoam | Med | Sintered Zirconia & Ceramic (Suitable for NP & TITN) | 0.12 | 6.5 | 10k/20k | \$28.75 Ea. |
| MCE210104C | Whl Flat Edge | Yellow | Crs | Sintered Zirconia & Ceramic (Suitable for NP & TITN) | 22 | 4.5 | 10k/20k | \$46.00 Ea. |
| MCE239104C | Disc | Yellow | Crs | Sintered Zirconia & Ceramic (Suitable for NP & TITN) | 22 | 2 | 10k/20k | \$37.00 Ea. |
| MCE239900M | Disc | Seafoam | Med | Sintered Zirconia & Ceramic (Suitable for NP & TITN) | 22 | 2 | 10k/20k | \$32.50 Ea. |
| CCE235900M | Lense | White | Med | Ceramic | 22 | 3.5 | 6K/20k | \$1.40 Ea. |

Diamond impregnated rubber and stones for zirconia and ceramics.

HANDPIECE POLISHERS AND ABRASIVES

For low speed handpiece systems




PSC138104F



PSC151104F



PSC199104F



PSC508104F



PSCSET104/01




TCP100104MSW4Z



TCP140104M




TCP158104MSW4Z



TCP546104M



TCP546104XF



TCPSET104/01



TCPSET104/02




PANTHER PREMIUM POLISHERS AND ABRASIVES

| Item No. | Shape | Color | Grit | Material Usage | Ø (mm) | [H] (mm) | Rpm Opt/Max | Price |
|----------------|----------------------|-----------------|------|------------------------------|--------|----------|-------------|-------------|
| PSC138104F | Cylinder | Panther Gray | Fine | Sintered Zirconia & Ceramic | 4 | 10.5 | 15k/20k | \$20.75 Ea |
| PSC151104F | Wheel Flat Edge | Panther Gray | Fine | Sintered Zirconia & Ceramic | 12 | 2.4 | 10k/20k | \$32.50 Ea |
| PSC199104F | Flat Cup | Panther Gray | Fine | Sintered Zirconia & Ceramic | 5 | - | 15k/20k | \$22.99 Ea |
| PSC508104F | Point (Rounded) | Panther Gray | Fine | Sintered Zirconia & Ceramic | 2.5 | 7 | 15k/20k | \$19.99 Ea |
| PSCSET104/01 | Panther Stone Kit | Panther Gray | Var | Sintered Zirconia & Ceramic | - | - | 15k/20k | \$142.99 Ea |
| TCP100104MSW4Z | Flame | Panther Pink | Med | Sintered Zirconia & Ceramic | 5.5 | 18 | 12k/20k | \$27.50 3Pk |
| TCP140104M | Disc | Panther Pink | Med | Sintered Zirconia & Ceramic | 26 | 2 | 8k/20k | \$44.99 Ea |
| TCP158104MSW4Z | Wheel Flat Edge | Panther Pink | Med | Sintered Zirconia & Ceramic | 17 | 2.4 | 8k/20k | \$28.00 Ea |
| TCP546104M | Spiral Polisher Disc | Panther Pink | Med | Spiral Polisher for Zirconia | .12 | 7 | | \$23.00 Ea |
| TCP546104XF | Spiral Polisher Disc | Panther White | XtFn | Spiral Polisher for Zirconia | .12 | 7 | | \$23.00 Ea |
| TCPSET104/01 | Panther Polisher Kit | Panther Wt/Pnk | Var | Sintered Zirconia & Ceramic | - | - | | \$288.99 Ea |
| TCPSET104/02 | Panther Abutment Kit | Panther Pnk/Yel | Var | Abutments | - | - | 8k/20k | \$157.00 Ea |

Panther precision stones and polishers provide smooth performance & incredible durability!

SUCTION SOLUTIONS

For individual workstations and milling machines



FZ VARIO CENTRAL SUCTION SYSTEMS

The premier solution for any laboratory. Powerful, clean and efficient; the advantages of central suction become obvious when the debris produced by the work stations or devices in the laboratory is reduced to an absolute minimum. Vacuum motors and filtration systems are available in different output classes and combinations allowing for a custom installation designed specifically for your laboratory. Suitable for providing suction at all dedicated workstations, for hand piece or resin processing, or for any dust producing device in dental laboratories of any size. Even devices with internal blower motors can be retrofit for connectivity to the central suction system.

For more information or a formal quotation, please contact: sales@zublerusa.com

FZ VARIO Central Suction System Parts

| Description | Item No. | Price |
|--|----------|----------|
| Liner Bags - FZ5-2,FZ10,FZ20,EF5, EF10 (Pkg. 50) | 995100 | \$102.00 |
| Filter Cartridge - FZ10,FZ20,EF10 | 556/052 | \$348.00 |
| Filter Cartridge - FZ5, FZ5-2. EFS | 556/062 | \$220.00 |
| Filter Cartridge - FZ5, FZ5-2. EFS PTFE | 556/062P | \$285.00 |



ZUBLER SINGLE & MULTI STATION SUCTION

Advanced Suction Systems



Z4 VARIOmatic
4 Workstations



Z1-CAM
For the milling machine



Z1-ECO
For high volume stations



Z1-M51

Z1-M50
For Single Station and some equipment



Z1-M11

Z1-M10

ZUBLER Z1 SUCTION UNITS

For Single Work Station, Equipment or Milling Machines

| Item# | Description | Price |
|----------|---------------------------------|------------|
| 800/0009 | Z1-AT Drawer Suction 120V | \$1,650.00 |
| 821/049 | Z1-CAM PRO Suction Unit 120V | \$3,659.99 |
| 821/048 | Z1-CAM Suction Unit 120V | \$2,397.00 |
| 821/047 | Z1-ECO PRO Suction Unit 120V | \$3,659.99 |
| 821/046 | Z1-ECO Suction Unit 120V | \$2,397.99 |
| 825/0136 | Z1-M10 Mobile Suction Unit 120V | \$1,939.99 |
| 825/0137 | Z1-M11 Suction Unit 120V | \$1,939.99 |

Quiet motors with legendary Zubler reliability and great suction capabilities.

ZUBLER MULTI STATION SUCTION SYSTEMS

For Multiple Work Stations, Equipment or Milling Machines


| Item# | Description | Price |
|----------|--|------------|
| 822/042 | Z4 VARIOmatic Austerity Package 4 Stations | \$8,590.00 |
| 822/0600 | AP600 | \$630.00 |
| 822/0605 | AP600-M | \$450.00 |
| 822/0615 | AP600-R | \$695.00 |
| 822/0612 | AP601 for Z4 incl. accessories | \$569.99 |

Use the Zubler app to take full advantage of all the features your suction system has to offer!



ZUBLER SUCTION ACCESSORIES


Parts and Accessories for Zubler Suction



ZUBLER Suction System Accessories

| Description | Item No. | Price |
|---|----------|----------|
| R1000 Rectangular pipe (Fit: R1200 & R1250) | 570/065 | \$83.99 |
| R1300 Retractable extension pipe (Fit: R1200 & R1250) | 825/228 | \$131.00 |
| R1200 Large Suction Funnel w/glass shield | 825/241 | \$204.99 |
| R1250 Small Suction Funnel w/glass shield | 825/236 | \$237.00 |
| R1500 (Flexible Suction Hood) | 825/273 | \$586.99 |

| Description | Item No. | Price |
|---|----------|----------|
| SH-1 Suction Hood System | 802/0035 | \$187.00 |
| Safety Glass Shield 140x190mm (Fit: R1200 & SH-1) | 557/172 | \$37.00 |
| Screen for SH-1 | 556/023 | \$22.00 |
| Monomer Suction Bath R1601 Plug-In | 825/229 | \$360.00 |
| RS 1 Cleaning Set | 825/490 | \$100.99 |
| RS 2 Cleaning Set | 825/497 | \$247.00 |





ZUBLER Suction Equipment Parts and Accessories


| Pic. | Description | Item No. | Price |
|------|--|-----------|----------|
| 1 | Carbon Brushes *New Style* | 037/10220 | \$21.00 |
| 2 | Carbon Brushes *Old Style* | 037/1010 | \$17.00 |
| 3 | Filter Bag - V6000 / V7000 | 556/035 | \$47.00 |
| 4 | Filter Bag - Z1 Eco / Z1 Cam | 556/040 | \$47.00 |
| 5 | Filter Bags (3 Pack) - V4010 / V4011 / Z1-M10 / Z1-M11 | 556/0035 | \$47.00 |
| 6 | FZ Filter Cartridge FZ10/20 | 556/052 | \$337.00 |
| 7 | FZ Liners Bags (Pkg. 50) | 995100 | \$102.00 |
| 8 | Micro Fine Filter (Cylinder) V Series / Z Series | 556/019 | \$69.00 |
| 9 | V-Series Vac Motor *New Style* | 450/008U | \$241.00 |

| Pic. | Description | Item No. | Price |
|------|--------------------------------|-----------|----------|
| - | Z1 Adapter - Exhaust Connector | 896/1103 | \$47.00 |
| 2 | Z1 Carbon Brushes Pair | 896/1102 | \$28.00 |
| - | Z1 Cable Adapter DC1-DC7 | 825/25617 | \$66.00 |
| - | Z1 Cord for Imes 3 Meter | 825/25616 | \$66.00 |
| - | Z1 Cord for Roland 3 Meter | 825/25613 | \$66.00 |
| - | Z1 Cord Wieland/VHF 3 Meter | 825/25610 | \$66.00 |
| - | Z1 Motor Only 100-120 V | 896/1105 | \$265.00 |
| - | Z1 PRO Filter Permanent | 556/0054 | \$102.00 |
| 10 | V-Series Vac Motor *Old Style* | 450/011 | \$241.00 |


OLD Style - For Serial Numbers Prior to: 21-7486-00

ZUBLER SUCTION ACCESSORIES

Extraction solutions for the entire laboratory



AV1000 WORKSTATION - A laboratory that strives for economic success in all its operations must keep a close eye on both quality and quantity. Usually this begins right in the preparatory stage for each process, where the foundation is laid for any top-quality work. This is why the design of the Zubler dust preparation box emphasizes ergonomics, cleanliness and efficiency.



AV1000 WORKSTATION

Preparation is the foundation for top-quality work. This is why the design of the AV1000 dust preparation box! It emphasizes ergonomics, cleanliness and efficiency.



KOTA Lighting

| DESCRIPTION | ITEM# | PRICE |
|---|-----------|----------|
| Lumild Luminaire | 424114093 | \$265.00 |
| Brighten your bench with this adjustable color correcting illumination. | | |

ZUBLER Stations & Accessories

| DESCRIPTION | ITEM# | PRICE |
|----------------------------|----------|------------|
| ZUBLER AV1000 Workstation | 825/310U | \$1,975.00 |
| FG 1000 Stand for AV1000 | 825/324 | \$495.00 |
| ZUBLER MB 1000 Monomer Box | 825/315 | \$1,395.00 |
| ZUBLER R1700 Spray Box | 825/223 | \$499.00 |



WASSERMANN Suction Box

| DESCRIPTION | ITEM# | PRICE |
|---------------------------|--------|------------|
| Suction Box Compact III | 143390 | \$1,145.00 |
| Arm Rests (Set) SBCIII | 143303 | \$Call |
| Cover Plates (Set) SBCIII | 143305 | \$Call |

WASSERMANN BOIL OUT UNITS

For Polymerization of multiple flasks



WAPO-EX12 II

WASSERMANN Boilout

| Description | Item No. | Price |
|---|----------|----------|
| WAPO-EX12 II | 170971 | \$Call |
| Boil out unit from Wassermann is a free standing unit for scalding up to 12 flask halves. | | |
| WAPO-EX 8 Plus | 170981 | \$Call |
| Compact tabletop boil out unit from Wassermann for scalding up to 8 flask halves. | | |
| WAPO-MAT III | 170980 | \$Call |
| Tabletop boil out unit from Wassermann for scalding up to 9 flasks in a holder at one time. | | |
| Wax Remover 2 Liter | 655001 | \$177.99 |

WAPO-EX 8 Plus




WAPO-Mat III



BDT EQUIPMENT

Bench & Removable Lab Equipment



BDT Equipment

| Description | Item No. | Price |
|--|----------|------------|
| Polymerization Bath IV0-4 | 910003 | \$3,595.00 |
| Plaster Separator | 910001 | \$319.00 |
| Plaster Separator Liner Bags Pkg of 25 | 920004 | \$28.00 |
| Extractor w/activated charcoal filter | 920192 | \$1,295.00 |
| Air Hose Return Device | 825/2642 | \$187.00 |
| Holder - straight with sub-structure | 825/2647 | \$21.00 |
| Holder - slanted with sub-structure | 825/2648 | \$21.00 |





Duplicating & Putty Silicone

| Description | Item No. | Price |
|---|-------------------|-------------------|
| Replisil (2 x 930mL) | 241122701 | \$99.00 |
| Addition-curing (1:1), Shore Hardness A 22 Curing-time 30 min | | |
| Stencil Putty Silicone (2 x 960mL) | 0261190215 | \$105.00 |
| Stencil Putty Silicone (5kg) | 0261190350 | \$290.00 |
| Stencil Putty Silicone (2 x 960mL) | 0261190425 | \$1,355.00 |

Lab Putty Addition-curing (1 :1) - Shore Hardness A 90

Cutting Discs, Glass Beads & Alumina Oxide

| Description | Item No. | Price |
|--|-----------|---------|
| Fiber Disc S (22mm x 0.3mm) (Pkg of 50) | 103220350 | \$95.00 |
| "S" = Small for all Crown & Bridge Non-Precious Alloys | | |
| White Alumina Oxide (50 micron/ 25lbs) | 995240W | \$72.99 |
| Glass Beads (50 micron / 25lbs) | 995413 | \$35.00 |


High quality abrasives last longer and clean faster.

MSDS and reference material available online at: www.zublerusa.com





sales@zublerusa.com



Atlanta, GA - (770) 921-2131
805 Marathon Parkway
Suite 130
Lawrenceville, GA 30046
8am- 5pm EST. Monday- Friday

Dallas, TX - (972) 600-9187
16135 Preston Rd.
Suite 126
Dallas, TX 75248
8am- 5pm CST. Monday- Friday