



Equipment and Materials for the Dental Professional



Over 65 years building reasons to
smile in more than 60 countries
worldwide

Product Guide
2023

NEW! PORTUX 3D PRINT RESIN



portux^{3D}

3D printing resin



PORTUX 3D PRINTING RESIN

Item #	Description	Price
P 111191	SHADE ORANGE x 1lb	\$30.00
P 0020104	SHADE BLACK x 1lb	\$30.00
P 111193	SHADE BLUE x 1lb	\$30.00
P 23972	SHADE FUCHSIA x 1lb	\$30.00
P 23971	SHADE GREEN x 1lb	\$30.00
P 23970	SHADE PINK x 1lb	\$30.00
P 0020103	SHADE WHITE x 1lb	\$30.00
P 23969	SHADE YELLOW x 1lb	\$30.00

COMPATIBLE PRINTERS	nm	PORTUX 3D Model	PORTUX 3D Model Ortho
ASIGA Pico 2	405	✓	✓
ASIGA MAX UV	385	✓	✓
Uniz Slash 2	405	✘	✘
Uniz IBEE	405	✓	✓
WhipMix VeriBUILD 3D	405	✓	✓
Flashforge Hunter	405	✘	✘
Anycubic Mono X	405	✓	✓
Anycubic Mono SE	405	✘	✓
Shining 3D Accufab-D1	405	✓	✓
Shining 3D Accufab-L4	405	✓	✓
Ackuretta Sol	405		✓
Phrozen Sonic Mini 4K	405	✓	✓
Phrozen Sonic XL 4K	405		✓
Rapidshape P30	385	✓	✘
Rapidshape P40	385	✓	✘

ORTHADONTIC ACRYLIC RESIN



O-Cryl Orthodontic Resin Monomer

O-cryl® is used for making brightly colored orthodontic appliances with exceptional adaptability and durability. O-cryl is highly resistant to fractures and allows an easy trimming and polishing to recover and maintain its original gloss.

OCRYL ORTHODONTIC RESIN

Item #	Description	Price	Item#	Description	Price	Item #	Description
P 111191	SHADE ORANGE x 1lb	\$30.00	P 23194	SHADE CLEAR x 1lb	\$30.00	P 23847	O-CRYL ORTHODONTIC RESIN KIT 8 X 40g + 2 X 55ml
P 0020104	SHADE BLACK x 1lb	\$30.00	P 24452	SHADE CLEAR x 5lb	\$88.00	Price:	\$66.50
P 111193	SHADE BLUE x 1lb	\$30.00	P 23195	SHADE CLEAR x 22lb	\$297.00		
P 23972	SHADE FUCHSIA x 1lb	\$30.00	M 17918	LIQUID SHADE BLUE X 8oz	\$19.00		
P 23971	SHADE GREEN x 1lb	\$30.00	M 17917	LIQUID SHADE YELLOW X 8oz	\$19.00		
P 23970	SHADE PINK x 1lb	\$30.00	M 17919	LIQUID SHADE RED X 8oz	\$19.00		
P 0020103	SHADE WHITE x 1lb	\$30.00	M 23196	LIQUID SHADE CLEAR X 8oz	\$19.00		
P 23969	SHADE YELLOW x 1lb	\$30.00	M 24453	LIQUID SHADE CLEAR X 32oz	\$40.00		



HIGH IMPACT & HEAT CURING ACRYLICS



Opti-Cryl High Impact Denture Base Polymer

- Strength to resist breakage
- Excellent handling characteristics
- Provides an accurate copy of the wax pattern
- Available in Original and Meharry shades
- Biocompatibility with oral tissues

Opti-Cryl Heat Curing Denture Base Polymer

- Excellent handling characteristics
- Biocompatibility with oral tissues
- Provides an accurate copy of the wax pattern
- 10 minutes working time
- Exceptional flexural strength

OPTI-CRYL® HIGH IMPACT DENTURE BASE POLYMER

Item #	Description	Price
P 108694	ORIGINAL X 1 lb	\$43.00
P 83268	ORIGINAL X 5 lb	\$205.50
P 83257	ORIGINAL X 22 lb	\$513.00
P 19537	MEHARRY X 1 lb	\$43.00
P 19517	MEHARRY X 5 lb	\$202.50
P 19540	MEHARRY X 22 lb	\$513.00

OPTI-CRYL® HEAT CURING DENTURE BASE POLYMER.

Item #	Description	Price
P 83272	ORIGINAL X 1 lb	\$25.00
P 83271	ORIGINAL X 5 lb	\$83.00
P 83238	ORIGINAL X 22 lb	\$297.00

OPTI-CRYL® HEAT CURING DENTURE BASE POLYMER

Item #	Description	Price
P 81152	MEHARRY X 1 lb	\$25.00
P 81153	MEHARRY X 5 lb	\$83.00
P 81196	MEHARRY X 22 lb	\$297.00
P 83760	CLEAR X 1 lb	\$25.00
P 81156	CLEAR X 5 lb	\$83.00
P 81199	CLEAR X 22 lb	\$298.00
P 82990	LIGHT PINK VEINED X 1 lb	\$25.00
P 82491	LIGHT PINK VEINED X 5 lb	\$83.00
P 81991	LIGHT PINK VEINED X 22 lb	\$297.00
P 82970	DARK PINK VEINED X 1 lb	\$25.00
P 81742	DARK PINK VEINED X 5 lb	\$83.00

HIGH IMPACT & HEAT CURING ACRYLICS



Opti-Cryl Self Curing Repair Polymer
Primarily for the repair of removable prostheses and the production of appliances, repairs and relines.

- Easy handling
- High gloss finish
- 4-6 minutes working time
- 10 minutes polymerization time at room temp.
- Optimal dimensional stability after repairs

Opti-Cryl Pour Acrylic Polymer
Used for manufacturing removable prostheses and also producing rebases and relines.

- Excellent handling characteristics
- Biocompatibility with oral tissues
- Provides an accurate copy of the wax pattern
- 1-minute mixing time
- Optimizes working time

OPTI-CRYL SELF CURING ACRYLIC (REPAIR POLYMER)

Item #	Description	Price
P 83269	ORIGINAL X 1 lb	\$27.00
P 83255	ORIGINAL X 5 lb	\$88.00
P 83250	ORIGINAL X 22 lb	\$297.00
P 81154	MEHARRY X 1 lb	\$27.00
P 81155	MEHARRY X 5 lb	\$88.00
P 24184	MEHARRY X 22 lb	\$297.00
P 81757	LIGHT PINK VEINED X 1 lb	\$27.00
P 81151	LIGHT PINK VEINED X 5 lb	\$88.00
P 81192	LIGHT PINK VEINED X 22 lb	\$297.00
P 110735	CLEAR X 1 lb	\$27.00
P 81157	CLEAR X 5 lb	\$88.00
P 81200	CLEAR X 22 lb	\$297.00
P 81746	DARK PINK VEINED X 1 lb	\$27.00

OPTI-CRYL POUR ACRYLIC POLYMER

Item #	Description	Price
P 22600	ORIGINAL X 1 LB	\$27.00
P 22602	ORIGINAL X 5 LB	\$\$\$8.00
P 22604	ORIGINAL X 22 LB	\$297.00
P 24182	MEHARRY X 1 LB	\$27.50
P 24183	MEHARRY X 5 LB	\$\$\$8.00
P 24184	MEHARRY X 22 LB	\$297.00
P 81744	LIGHT PINK VEINED X 1 lb	\$27.50
P 81743	LIGHT PINK VEINED X 5 lb	\$\$\$8.00
P 81193	LIGHT PINK VEINED X 22 lb	\$297.00





OPTI-CRYL HEAT CURING ACRYLIC LIQUID MONOMER

Item #	Description	Price
80155	OPTI-CRYL HEAT CURING LIQUID X 8 oz	\$16.75
80156	OPTI-CRYL HEAT CURING LIQUID X 32 oz	\$40.00
80020	OPTI-CRYL SELF CURING LIQUID X 8 oz	\$16.75
80021	OPTI-CRYL SELF CURING LIQUID X 32 oz	\$40.00
108745	OPTI-CRYL POUR LIQUID X 8 oz	\$16.75
81198	OPTI-CRYL POUR LIQUID X 32 oz	\$40.00

OPTI-CRYL HEAT CURING ACRYLIC CHARACTERIZATION KIT contains: 8 x 40g bottles.

Colors: Flat Pink, Smooth C, Smooth Original, Smooth B, Smooth VR5, Smooth Meharry, Red V51 and Beige -A4.

OPTI-CRYL HEAT CURING ACRYLIC CHARACTERIZATION KIT

Item #	Description	Price
P 21953	OPTI-CRYL HEAT CURING CHARACTERIZATION KIT 8 X 40g	\$83.00
P 23428	OPTI-CRYL FIBERS FOR CHARACTERIZATION X 2 (Blue and Red)	\$12.00
P25708	OPTI-CRYL RED V51-HEAT CURE 125gm	\$30.00
P25710	OPTI-CRYL SMOOTH MEHARRY – HEAT CURE 125gm	\$30.00
P25712	OPTI-CRYL SMOOTH VR5 – HEAT CURE 125gm	\$30.00
P25714	OPTI-CRYL SMOOTH B – HEAT CURE 125gm	\$30.00
P25716	OPTI-CRYL SMOOTH ORIGINAL – HEAT CURE 125gm	\$30.00
P25718	OPTI-CRYL SMOOTH C – HEAT CURE 125gm	\$30.00
P25720	OPTI-CRYL FLAT PINK X – HEAT CURE 125gm	\$30.00
P25722	OPTI-CRYL BEIGE A4 – HEAT CURE 125gm	\$30.00

CROWN & BRIDGE TOOTH SHADE ACRYLIC RESINS

Novacryl L has all of the characteristics of a premium laboratory acrylic resin for temporaries. Novacryl L is a co-polymer allowing significantly more color stable for temporary crowns, bridges, inlays, and onlays. Ideal for laboratory use.



NOVACRYL SELF CURING ACRYLIC

Item #	Shade x Size	Price
P 17902	A0 X 125g	\$28.00
P 17890	A1 X 125g	\$28.00
P 17891	A2 X 125g	\$28.00
P 17903	A3 X 125g	\$28.00
P 17892	A3.5 X 125g	\$28.00
P 17893	A4 X 125g	\$28.00
P 17906	B0 X 125g	\$28.00
P 17904	B1 X 125g	\$28.00
P 17894	B2 X 125g	\$28.00
P 17895	B3 X 125g	\$28.00
P 17896	B4 X 125g	\$28.00

NOVACRYL SELF CURING ACRYLIC CONT.


Item#	Shade x Size	Price
P 17897	C1 X 125g	\$28.00
P 17900	C2 X 125g	\$28.00
P 17898	C3 X 125g	\$28.00
P 17905	C4 X 125g	\$28.00
P 17899	D2 X 125g	\$28.00
P 17901	D3 X 125g	\$28.00
P 17889	D4 X 125g	\$28.00
P 17908	INCISAL DARK X 125g	\$28.00
P 17907	INCISAL LIGHT X 125g	\$28.00
NOVACRYL LIQUID		
M 25707	LIQUID 8 OZ	\$16.75

PORTUX PMMA CAD/CAM DISCS



PORTUX Cad/Cam PMMA DENTURE & CASTABLE DISCS

Item #	Description	Price
NW22916	Clear Appliances/Castable 98.5 x 20	\$33.50
NW22917	Red Castable 98.5 x 20	\$33.50
NW22160	Denture Base Original	\$44.00
NW26110	Portux Denture Base Shade Meharry 98.5 x 30	\$44.00



PORTUX Cad/Cam PMMA MONO SHADED DISCS

Item #	SHD	SIZE	Retail	Item #	SHD	SIZE	Retail	Item #	SHD	SIZE	Retail	Item #	SHD	SIZE	Retail
NW21936	A1	98.5 X 18	\$36.50	NW21945	BL2	98.5 X 18	\$36.50	NW24139	B1	98.5 X 30	\$41.50	NW23324	C2	98.5 X 18	\$36.50
NW21924	A1	98.5 X 20	\$38.50	NW21933	BL2	98.5 X 20	\$38.50	NW24145	B1	98.5 X 35	\$46.50	NW23325	C2	98.5 X 20	\$38.50
NW22317	A1	98.5 X 25	\$40.25	NW22332	BL2	98.5 X 25	\$40.25	NW22371	B2	98.5 X 18	\$36.50	NW23326	C2	98.5 X 25	\$40.25
NW24136	A1	98.5 X 30	\$41.50	NW24141	BL2	98.5 X 30	\$41.50	NW22376	B2	98.5 X 20	\$38.50	NW24930	C2	98.5 X 30	\$41.50
NW24142	A1	98.5 X 35	\$46.50	NW24147	BL2	98.5 X 35	\$46.50	NW22381	B2	98.5 X 25	\$40.25	NW23333	D2	98.5 X 18	\$36.50
NW21939	A2	98.5 X 18	\$36.50	NW26536	BL3	98.5 X 18	\$36.50	NW24140	B2	98.5 X 30	\$41.50	NW23334	D2	98.5 X 20	\$38.50
NW21927	A2	98.5 X 20	\$38.50	NW26537	BL3	98.5 X 20	\$38.50	NW23321	C1	98.5 X 18	\$36.50	NW23335	D2	98.5 X 25	\$40.25
NW22320	A2	98.5 X 25	\$40.25	NW26538	BL3	98.5 X 25	\$40.25	NW23322	C1	98.5 X 20	\$38.50	NW24940	D2	98.5 X 30	\$41.50
NW24137	A2	98.5 X 30	\$41.50	NW26539	BL3	98.5 X 30	\$41.50	NW23323	C1	98.5 X 25	\$40.25	NW23337	D3	98.5 X 20	\$38.50
NW24143	A2	98.5 X 35	\$46.50	NW26540	BL3	98.5 X 35	\$46.50	NW24938	C1	98.5 X 30	\$41.50	NW23338	D3	98.5 X 25	\$40.25
NW22093	A3	98.5 X 18	\$36.50	NW21942	B1	98.5 X 18	\$36.50	<p style="text-align: center;">NEW!</p> <p style="text-align: center;">BL3 NOW AVAILABLE IN 18mm, 20mm, 25mm, 30mm & 35mm!</p>							
NW22088	A3	98.5 X 20	\$38.50	NW21930	B1	98.5 X 20	\$38.50								
NW22323	A3	98.5 X 25	\$40.25	NW22329	B1	98.5 X 25	\$40.25								
NW24138	A3	98.5 X 30	\$41.50												



Portux® multilayer discs are designed with a strategic distribution of 10 layers, providing incomparable natural-looking aesthetics resulting from the seamless transition of shade - cervical to incisal with natural translucency.

PORTUX Cad/Cam PMMA MULTILAYER DISCS

Item #	SHD	SIZE	Retail	Item #	SHD	SIZE	Retail	Item #	SHD	SIZE	Retail	Item #	SHD	SIZE	Retail
NW21846	A1	98.5 X 18	\$71.50	NW23118	A3	98.5 X 30	\$86.00	NW21852	B1	98.5 X 18	\$71.50	NW24477	C1	98.5 X 20	\$75.00
NW21858	A1	98.5 X 20	\$75.00	NW22218	A3.5	98.5 X 20	\$75.00	NW21864	B1	98.5 X 20	\$75.00	NW24478	C1	98.5 X 25	\$82.00
NW22145	A1	98.5 X 25	\$82.00	NW22225	A3.5	98.5 X 25	\$82.00	NW22154	B1	98.5 X 25	\$82.00	NW24868	C1	98.5 X 30	\$86.00
NW23116	A1	98.5 X 30	\$86.00	NW21855	BL2	98.5 X 18	\$71.50	NW23119	B1	98.5 X 30	\$86.00	NW24479	C2	98.5 X 18	\$71.50
NW22576	A1	98.5 X 35	\$92.00	NW21867	BL2	98.5 X 20	\$75.00	NW22577	B1	98.5 X 35	\$92.00	NW24480	C2	98.5 X 20	\$75.00
NW21849	A2	98.5 X 18	\$71.50	NW22157	BL2	98.5 X 25	\$82.00	NW22350	B2	98.5 X 18	\$71.50	NW24481	C2	98.5 X 25	\$82.00
NW21861	A2	98.5 X 20	\$75.00	NW23121	BL2	98.5 X 30	\$86.00	NW22355	B2	98.5 X 20	\$75.00	NW24870	C2	98.5 X 30	\$86.00
NW22148	A2	98.5 X 25	\$82.00	NW22579	BL2	98.5 X 35	\$92.00	NW22360	B2	98.5 X 25	\$82.00	NW24482	D2	98.5 X 18	\$71.50
NW23117	A2	98.5 X 30	\$86.00	NW27144	BL3	98.5 X 18	\$71.50	NW23120	B2	98.5 X 30	\$86.00	NW24483	D2	98.5 X 20	\$75.00
NW24457	A2	98.5 X 35	\$92.00	NW27145	BL3	98.5 X 20	\$75.00	NW24476	C1	98.5 X 18	\$71.50	NW24484	D2	98.5 X 25	\$82.00
NW22083	A3	98.5 X 18	\$71.50	NW27146	BL3	98.5 X 25	\$82.00					NW24872	D2	98.5 X 30	\$86.00
NW22078	A3	98.5 X 20	\$75.00	NW27147	BL3	98.5 X 30	\$86.00								
NW22151	A3	98.5 X 25	\$82.00	NW27148	BL3	98.5 X 35	\$92.00								

NEW! BL3 Now available in 18mm, 20mm, 25mm, 30mm & 35mm!



REFLECTION

Acrylic teeth are designed to replace the loss of one or several dental pieces. They are made of Polymethyl Methacrylate, Pigments, Ethylene glycol dimethacrylate -cross Link-, Fluorescence, and Methyl Methacrylate. Acrylic teeth come in a wide variety of tooth shades and forms that can be adapted to each patient's needs. New Stetic has different marks for its acrylic teeth according to their number of layered structures.

The gingival layer is the one that gives the tooth color. This layer is placed inside the tooth and simulates the tooth's dentine beneath the enamel; the incisal layer is the external tooth's layer that simulates the enamel. This layer gives translucence and veined effect to artificial teeth.

FOUR-LAYER ACRYLIC TEETH: Two gingival layers + two incisal layers


Anteriors teeth: upper 9 and lower 5

Posteriors teeth: upper and lower 4

Shades: 18

A0, A1, A2, A3, A3.5, A4, B0, B1, B2, B3, B4, C1, C2, C3, C4, D2, D3, D4

MICRO-FLASK was developed to maximize the results in microwave curing procedures. Made of a Reinforced Engineered Plastic Formula that guarantees the best physical and mechanical properties in the market. Our Flask is specifically designed to perfection to be used with OPTI-CRYL POUR as well with other brands in the market.



NEWSTETIC ACCESSORIES

Item #	Description	Price
F 0016779	MICROWAVE POLYMERIZATION FLASK	\$250.00
C 17934	#2 MEDIUM COTTON ROLL	\$36.00
K 97190	KROMALGIN BAG X 453 g	\$19.50
C 0015336	PLASTIC CABINET X 10 DRAWERS	\$166.50
0025700	PRINTED SHADE GUIDE	\$0.00
G 111271	O-CRYL SELF CURING ACRYLIC SHADE GUIDE X 8+3 TABLETS	\$0.00
G 93413	OPTI-CRYL ACRYLIC SHADE GUIDE X 6 TABLETS	\$0.00
P 15822	OPTI-CRYL SELF CURING ACRYLIC TRAYS BLUE X 1lb	\$16.50
D 11188	DURACRYL FLOW INLAY RESIN FOR PATTERNS 60g+55cc	\$35.00

zubler USA

 ceramay[®]
dermalceramics

 mihmvoigt

 zubler[®]

 APPLIEDENTAL

 WASSERMANN
Hamburg · Germany

 panther

 newstetic
USA

 BDT

 DENTAL
Concept Systems

nssales@zublerusa.com



Atlanta, GA - (770) 921-2131
805 Marathon Parkway
Suite 130
Lawrenceville, GA 30046
8am- 5pm EST. Monday- Friday

Dallas, TX - (972) 600-9187
16135 Preston Rd.
Suite 126
Dallas, TX 75248
8am- 5pm CST. Monday- Friday