




VARIOstar T50

Tabletop handpiece system



Distributed by:

zubler USA

Zubler USA
130 Satellite Blvd NE, Ste E
Suwanee, GA 30024

Office (770) 921-2131
www.zublerusa.com





System-Variants:

804 - 3112 VARIOstar K50 Knee control unit 220 V - 240 V

804 - 3111 VARIOstar K50 Knee control unit 115 V

Contains:

804 - 3012 Knee control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3011 Knee control unit 115 V



804 - 3152 VARIOstar T50 Tabletop unit with Foot pedal 220 V - 240 V

804 - 3151 VARIOstar T50 Tabletop unit with Foot pedal 115 V

Contains:

804 - 3052 Tabletop control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3051 Tabletop control unit 115 V

804 - 1016 Foot pedal



804 - 3162 VARIOstar T50 Tabletop unit with Foot control lever left/right 220 V - 240 V

804 - 3161 VARIOstar T50 Tabletop unit with Foot control lever left/right 115 V

Contains:

804 - 3052 Tabletop control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3051 Tabletop control unit 115 V

804 - 1015 Foot control lever



804 - 3172 VARIOstar T50 Tabletop unit with Knee pedal 220 V - 240 V

804 - 3171 VARIOstar T50 Tabletop unit with Knee pedal 115 V

Contains:

804 - 3052 Tabletop control unit 220 V - 240 V

804 - 3002 Handpiece with chuck Ø 2,35 mm

804 - 3051 Tabletop control unit 115 V

804 - 1017 Knee Pedal (left/right mountable)



All systems include a tool for changing the collet chuck, a handpiece holder and operating instructions.

Technical data:

Length: 162 mm
Weight: 235 g
Rated voltage: AC 100V/115V/230V/240V
Rated frequency: ~ 50/60Hz
Speed: 1000-50000 rpm
Chuck: Ø 2.35 mm or optional Ø 3.00 mm
Torque: 7,8 Ncm
Power consumption: max. 270 W
Protection class: I

Variants: Knee control unit
Tabletop unit with knee pedal
Tabletop unit with foot pedal
Tabletop unit with foot control left/right
Warranty: Lifetime warranty on armature and winding
All other parts 1 Year

Subject to technical changes.

Distributed by:

zubler USA

Zubler USA
130 Satellite Blvd NE, Ste E
Suwanee, GA 30024

Office (770) 921-2131
www.zublerusa.com

