

**#1 SELLING HANDPIECE
IN THE WORLD!**

OVER 1 MILLION SOLD!



MODEL: A1

CONTENTS

01

Features

02

Functional Description

03

Advantages

04

Product Parameters



FEATURES

MODEL:A1

- Push button
- Torque handpiece
- Clean head system
- Four points anti-retraction
- Single water spray



FEATURES



Updated Peacock grip
design



Japanese Ceramic
Bearing



Patented Spindle
Design for secure bur
retention.

FEATURES

JAPANESE CERAMIC BEARINGS




www.zublerusa.com

FEATURES



1. Clean head system - To prevent cross-infection



2. Super torque - For efficient Drilling

sales@zublerusa.com

TECHNICAL DATA

Noise : $\leq 60\text{dB}$

Vibartion : $\leq 0.04\text{mm}$



TECHNICAL DATA

COMPARISON

	Appledental	Being	COXO
Noise:	≤60dB	≤65dB	≤65dB
Torque/Power:	22W-24W	20W-22W	20W-22W
Warranty:	6 months	6 months	6 months

sales@zublerusa.com

TECHNICAL DATA

Code:	A1
Air Pressure:	0.25 (Mpa) -0.3 (Mpa)
Rotor Speed:	320 (krpm) -400 (krpm)
Torque/Power:	22W-24W
Water flow:	90(min/ml)-110(min/ml)
Water pressure:	0.2 (Mpa)-0.23 (Mpa)
Atomization pressure:	0.3 (Mpa)-0.35 (Mpa)
Noise:	≤60dB
Description:	torque handpiece,push button ,clean head system,Four points anti-retraction ,single water spray,4 holes and 2 holes optional,exposed cartridge

sales@zublerusa.com



Distributed by:

 **zubler** USA

130 Satellite Blvd NE
Suite E
Suwanee, GA 30024
sales@zublerusa.com
(770) 921-2131

www.zublerusa.com